

Digital SLR



Digital SLR Camera, 18 megapixel, Canon. **EOS Rebel T5i** \$779.00
 Includes Canon 18-55mm EF-S IS II lens, SD-16 GB card and shoulder bag.
AC Adapter Kit. ACK-E8 \$104.00


Accessories



Adapter, CS-Mount to C-Mount, 5 mm spacer \$15.00
Do not use this adapter ring on an existing C-mount camera, as it will negatively impact the quality and parfocality of the imaging system.

Cables

BNC, male-to-male, 75Ω, RG59/U
 6' \$8.80
 12' \$10.60
S-VHS, 4-pin male to 4-pin male
 6' \$10.00
 12' \$14.00




Cable Connectors
 BNC Female to RCA Male. **8110** \$5.65
 BNC Female to BNC Female. **25-7430** \$4.20
 BNC Tee, female-male-female. **27-8140** \$9.00
 RCA Female to BNC Male. **7510** \$4.50
Cable Release. AR-3 \$24.95



Couplers, Diagnostic Instruments, USA
 1X C-Mount for 38 mm ISO Photoport. **D10NLC** \$135.00
Fits newer Nikon and Leica microscopes
 Nikon SLR Lens to C-Mount TV. **NIKL-CMTC** \$116.00

Couplers, National
 1X C-Mount for 23.2 mm ID Eyetube. **930-023** \$29.00
Lens free - this coupler is not parfocal

Magnifiers



Loupes
 Eye, plastic, 10X at 1", black. **ELP-174.10** \$2.60
 Triplet, 10X at 1", 18 mm, black. **ELP-750.01** \$15.95
 Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.



Large diameter (4 1/8") glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

Magnifier, Folding (Linen Tester), 3X. ELP-800.03 \$9.95



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

Magnifier, Forceps. M15 \$4.00



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

Magnifier, Handle, Ultra-Optix
 2" Lens, 3X. **SV-2P** \$5.45
 3" Lens, 2X with 5X bifocal inset. **SV-3P** \$7.45
 4" Lens, 2X with 4X bifocal inset. **SV-4P** \$12.55



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

Magnifier, Headband, Optivisor®, USA

- 1.75X at 14", DA3. **ELP-503.00** \$48.95
- 2.5X at 8", DA5. **ELP-505.00** \$48.95
- 2.75X at 6", DA7. **ELP-507.00** \$50.95
- Light Attachment, Quasar™ LED. **ELP-558.00** \$29.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

Magnifier, Lighted, 2 x 2" lens, 2X. SV-BPL \$7.45



Confirm patient identification and preview the slide before examining microscopically. This macro inspection magnifier mounts onto the microscope's arm, where the lens angle can be adjusted using the double-ball joint arm. Large 3 x 1 3/4" lens, 2.25X, permits viewing of entire 3 x 1" glass slide.

Magnifier, Magnistitch™, USA. 15401 \$12.50