


**Rongeurs & Bone Instruments**




Blumenthal Rongeur

- Rongeurs, stainless steel**  
 Blumenthal, 15 cm (5¾")  
 Economy. **410-571** ..... \$41.30  
 Standard. **650-571** ..... \$112.00  
 Blumenthal, Slightly Angled Jaws, 15 cm (5¾"). **9571** ..... \$37.50  
 Mini-Blumenthal, 11.5 cm (4½"). **410-560** ..... \$53.10  
 Friedman, 14 cm (5½")  
 Economy. **410-574** ..... \$41.30  
 Standard. **650-574** ..... \$112.00




Liston Bone Cutter


- Bone Cutters, stainless steel, economy**  
 Liston  
 14 cm (5½"), 20 mm blades. **320-086** ..... \$24.30  
 21.5 cm (8½"), 30 mm blades. **320-087** (1095) ..... \$29.80  
 Scissor Style, 10 cm (4"), 18 mm blades. **320-025** ..... \$8.80  
 straight, strong tips




- Bone Nipper, 10 cm (4") long, 6 mm cutting edge. B286** ..... \$12.00  
 angled offset jaws



- Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws**  
 Economy. **33-1001** ..... \$45.00  
 German. **965-620** ..... \$160.00



- Postmortem Hammer, hook handle, 25 cm (10") long**  
 Economy. **33-1032** ..... \$35.70  
 Standard. **74-043** ..... \$80.00



16 mm (5/8") cutting blade;  
 9 cm (3½") long

- Skull Breaker, Virchow, economy. 33-1011** ..... \$42.40

**Scalpels**



Economy

**Scalpels**

- Disposable, plastic handle, sterile  
 Economy, gray handle
- |   | bx/10  | cs/50 bx |
|---|--------|----------|
| ★ #10 carbon steel blade. <b>9410</b> ..... | \$9.95 | 431.00   |
| #11 carbon steel blade. <b>9411</b> .....   | \$9.95 | 431.00   |
| #22 carbon steel blade. <b>9422</b> .....   | \$9.95 | 431.00   |
- Lance, green handle, England  
 #10 carbon steel blade, Bx/10. **92410** ..... \$13.00  
 #11 carbon steel blade, Bx/10. **92411** ..... \$13.00  
 #21 carbon steel blade, Bx/10. **92421** ..... \$13.80  
 #22 carbon steel blade, Bx/10. **92422** ..... \$13.80



- Miniature, plastic thumb handle, nonsterile, England  
 #11P stainless steel blade, each. **2201111** ..... \$0.95




Smooth Handle



Fluted Handle

- Reusable, one-piece stainless steel, 1½" blade  
 Smooth Handle, 15.5 cm (6"). **14-350** ..... \$3.15  
 (1445)  
 Fluted Handle, narrow blade, 15.5 cm (6"). **14-360** ..... \$3.15  
 (1460)

Continued



**Scalpel Blades**, individually wrapped

Carbon Steel, Economy, Bx/100 .....	\$19.65
#10 blade. <b>9310</b>	#21 blade. <b>9321</b>
#11 blade. <b>9311</b>	#22 blade. <b>9322</b>
#12 blade. <b>9312</b>	#23 blade. <b>9323</b>
#15 blade. <b>9315</b>	#24 blade. <b>9324</b>
#20 blade. <b>9320</b>	
Carbon Steel, Bard-Parker, Rib-Back® design, nonsterile, <b>USA</b>	
#60 blade, 25 strips of 6, Bx/150. <b>371340</b> .....	\$343.00
Carbon Steel, Germany, nonsterile	
#60 blade, 25 strips of 4, Bx/100. <b>2200260</b> .....	\$46.00
#70 blade, 25 strips of 4, Bx/100. <b>2200270</b> .....	\$63.00
Carbon Steel, Swann-Morton, nonsterile	
#60 blade, Bx/50. <b>22002SM60</b> .....	\$29.60
Stainless Steel, Economy, Bx/100 .....	\$21.95
#10 blade. <b>36-010S</b>	#21 blade. <b>36-021S</b>
#11 blade. <b>36-011S</b>	#22 blade. <b>36-022S</b>
#20 blade. <b>36-020S</b>	
Stainless Steel, Lance, Bx/100 .....	\$38.00
#11 blade. <b>92111</b>	#15 blade. <b>92115</b>
Stainless Steel, Economy	
#60 blade, Bx/100. <b>052-620</b> .....	\$45.00
#70 blade, Bx/100. <b>052-622</b> .....	\$45.00



Safely remove and store up to 100 used blades in this disposable box.

**Scalpel Blade Remover. 22048** .....

\$3.60



**Scalpel Handles**, reusable

Stainless Steel, solid	
#3 Style, accepts blades #10-17	
Economy. <b>36-405</b> .....	\$4.20
Standard. <b>10-305</b> .....	\$11.50
#4 Style, accepts blades #20-26	
Economy. <b>36-410</b> .....	\$4.20
Standard. <b>10-310</b> (565-410) .....	\$11.50
#7 Style, accepts blades #10-17	
Economy. <b>36-415</b> .....	\$4.20
Standard. <b>10-315</b> (565-415) .....	\$18.50

Continued

Stainless Steel, Hollow

#8 Style, accepts blades #60 & #70. <b>HAN8</b> .....	\$8.00
PVC Coated Stainless Steel. <i>Not autoclavable</i>	
#3 Style, accepts blades #10-17. <b>KNF-225.00</b> .....	\$5.35
Plastic, contoured to fit palm of hand	
#5 Style, accepts blades #10-17. <b>KNF-221.00</b> .....	\$5.00
#6 Style, accepts blades #20-26. <b>1605</b> .....	\$5.00
#8 Style, accepts blades #60 & #70	
Economy. <b>HAN8P</b> .....	\$6.00
Bard-Parker™, <b>USA. 371080</b> .....	\$31.35

**Scissors**




Barraquer


**Scissors, Spring**, ambidextrous, stainless steel

Barraquer (butterfly), 6 cm long, economy. <b>3737</b> .....	\$34.50
10 mm angled-on-flat blades	

36205/5-254/5-855



M3720



Castroviejo, straight, sharp tips, economy


7 cm (2¾"), 5 mm blades. <b>36205</b> .....	\$45.00
Effective cutting edge is 3 mm; tip-to-tip is 3 mm.	
8 cm (3"), 7 mm blades. <b>5-254</b> .....	\$45.00
Effective cutting edge is 6 mm; tip-to-tip is 4 mm.	
8.5 cm (3¼"), 5 mm blades. <b>5-855</b> .....	\$45.00
Effective cutting edge is 4 mm; tip-to-tip is 3 mm.	
11.5 cm (4½"), 17 mm blades. <b>M3720</b> .....	\$26.00
Effective cutting edge is 12 mm.	



Uniband™

Lifetime guarantee against rust and defects.

Featherlite™



Uniband™, sharp tips, Laschal Surgical, **USA**

Straight, 12 mm Blades. <b>LA-4</b> .....	\$32.00
Straight, 20 mm Blades. <b>LA-1</b> .....	\$32.00
Straight, 20 mm Serrated Blades. <b>LA-2</b> .....	\$32.00
Curved, 12 mm Blades. <b>LA-4C</b> .....	\$32.00
Curved, 20 mm Blades. <b>LA-1C</b> .....	\$32.00

Featherlite™, Castroviejo, sharp tips, Laschal Surgical, **USA**

Greater tactile sensitivity than Uniband™ scissors.

Straight, 12 mm Blades. <b>MPF-4</b> .....	\$80.80
Straight, 20 mm Blades. <b>MPF-1</b> .....	\$80.80
Curved, 12 mm blades. <b>MPF-4C</b> .....	\$80.80

**Scissors, Bonn**, 15 mm blades, 9 cm (3½"), Germany

Straight, Roboz. <b>RS-5840</b> .....	\$89.10
Curved, Roboz. <b>RS-5841</b> .....	\$90.10