

**Inverted Microscopes**



CP-2A1

Jenco

**CP-2A1** Inverted, Trinocular, Brightfield/Phase..... \$2,975.00

- Widefield 10X/20 eyepieces
- Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Six DIN, LWD plan achromat objectives: brightfield 10X, 25X, and 40X; plan phase 10X, 25X, and 40X
- 5-place, reverse facing nosepiece
- Attachable mechanical stage includes three inserts
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
- Centerable, dedicated phase annuli mounted on sliders
- 30 watt halogen Koehler-type illuminator
- 48.25 cm (19")D x 21.6 cm (8½")W x 47 cm (18½")H
- 5-year mechanical & optical/2-year electrical warranty



BC-366

Jenco

**BC-366** Inverted, Trinocular, BF/Phase/Fluor ..... \$8,960.00

- Widefield 10X/22 eyepieces
- Binocular, 45° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Side camera port
- Six infinity, LWD plan achromat objectives: brightfield/fluorescence 10X, 25X, and 40X; phase 10X, 25X, and 40X
- 6-place, angled nosepiece
- Attachable mechanical stage
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; centerable, dedicated phase annuli mounted in turret
- 30 watt halogen transmitted light; Koehler-type illuminator
- 100 watt mercury epi-fluorescence
- Blue (FITC) and green (TRITC) fluorescence filter sets
- 5-year mechanical & optical/2-year electrical warranty

Meiji ..... Access [www.meijitechno.com](http://www.meijitechno.com)  
and call us to place an order

**Polarized Light Microscopes**



Meiji ..... Access [www.meijitechno.com](http://www.meijitechno.com)  
and call us to place an order

**Metallurgical Microscopes**



Jenco ..... Access [www.jencointernational.com](http://www.jencointernational.com)  
and call us to place an order

Meiji ..... Access [www.meijitechno.com](http://www.meijitechno.com)  
and call us to place an order