


Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. 1596..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CBP12** \$13.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. DEM-345.00 \$80.95



Tools not included.

Dissecting Cases

6-Loop with snap fastener, black leatherette. **1257370**..... \$1.15
 12-Loop with snap fastener, black leatherette. **1257371**..... \$1.60
 13-Loop with zipper, brown leatherette. **36900** \$7.80
 Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits


*Inquire about customizing a kit for your surgical procedure.
 Kit components are discounted 10%.*

Basic kit. 38501..... \$5.30
 Contents: Dissecting Case, 6-loop (7-370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551..... \$17.85
 Contents: Dissecting Case, 6-loop (7-370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20 7/8 x 12 13/16 x 2 1/2". **AM20029** \$21.90
 25 3/4 x 17 3/4 x 2 1/4". **AM68357** \$64.90



Dissecting Pins, T-pin, 2", nickel-plated steel

1/2 lb Box, ~430 pcs. **214500**..... \$14.85



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic

Black Enameled, Pk/100 \$7.00
 #000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
 #00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
 #0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**
 Stainless Steel, Pk/100 \$8.75
 #000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
 #00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
 #0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. 150-1005 \$43.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1 1/2". **35-406** \$26.60
 Large, 10 x 4 x 2". **35-408** \$35.00



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450..... \$419.00
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455** \$21.00

Continued