




SLIDE & CASSETTE PROCESSING

Microscope Slides

Slides, general, 75 x 25 x 1 mm

	bx/72	cs/20 bx
Plain, ground edges		
Premier®, safety corners. 9101 	\$3.90	76.00
 Globe, 90° corners, "Globe". 1301	—	99.00
Frosted, beveled edges		
Premier®, sandblasted, safety corners. 8205	\$6.05	115.00
 Globe, etched, clipped corners, "Globe". 1334 ...	—	120.00
White End, clipped corners, "Premiere"		
Premier®. 9108W-E	—	115.00
Superfrost®, "SPECIMEN". 4951-001T	—	257.00
(12-550-143, 48311-600)		

Superfrost® Plus



Concavity, 2-Well



Slides, specialty

Charged (silanated), 75 x 25 x 1 mm
Recommended for adhering tissue sections
and cell preparations.


Diamond™, white end, Globe		
Cs/20 Bx. 1358W	\$264.00	
Superfrost® Plus, white end, Erie, USA		
Cs/20 Bx. 4951PLUS-001	\$597.00	
Preferred slide for Ventana instrumentation. (12-550-15, 48311-703)		


Concavity


75 x 25 x 1.2 mm, economy	ea	bx/12
1-Well, 15 mm dia., 0.5 mm deep. 7103	\$0.70	6.30
2-Well, 15 mm dia., 0.5 mm deep. 7104	\$0.95	8.45
76 x 26 x 3.2 mm, Erie		bx/12
1-Well, 18 mm dia., 1.75 mm deep. 1526	\$64.95	
2-Well, 18 mm dia., 1.75 mm deep. 1527	\$72.55	
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA		
QMS3X1X1 , each.....	\$20.00	
Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA		
FSMS3X1X1 , each.....	\$29.00	

Coverslips (Cover Glasses)


Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)		
Circles, Germany* 		
25 mm, ~222 pcs/oz, Cs/10 oz. 0-25		\$214.00



#1 Thickness (0.13-0.17 mm)		
Squares & Rectangles, Erie, USA 	1 oz	cs/10 oz
22 x 22 mm, ~162 pcs/oz. 1-2222	\$7.45	62.70
22 x 30 mm, ~119 pcs/oz. 1-2230	—	57.30
22 x 50 mm, ~71 pcs/oz. 1-2250	—	57.30
24 x 50 mm, ~65 pcs/oz. 1-2450	—	57.30
24 x 60 mm, ~54 pcs/oz. 1-2460	—	57.30

#1 Thickness (0.13-0.17 mm)		
Circles, Germany* 	cs/5 oz	cs/10 oz
12 mm, ~663 pcs/oz. 1-12	\$268.00	—
18 mm, ~296 pcs/oz. 1-18	\$150.00	—
25 mm, 100 pcs/pk, Cs/10 Pk. 1-25		\$212.00

#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*

Squares & Rectangles, Erie, USA 	1 oz	cs/10 oz
18 x 18 mm, ~209 pcs/oz. 1.5-1818	—	61.40
22 x 22 mm, ~140 pcs/oz. 1.5-2222	\$7.45	62.70
22 x 40 mm, ~77 pcs/oz. 1.5-2240	—	57.30
22 x 50 mm, ~62 pcs/oz. 1.5-2250	—	57.30
24 x 30 mm, ~94 pcs/oz. 1.5-2430	—	57.30
24 x 40 mm, ~71 pcs/oz. 1.5-2440	—	57.30
24 x 50 mm, ~56 pcs/oz. 1.5-2450	\$6.85	57.30
24 x 60 mm, ~47 pcs/oz. 1.5-2460	—	57.30

#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*

Circles, Germany* 	cs/5 oz	cs/10 oz
12 mm, ~552 pcs/oz. 1.5-12	\$223.00	—
18 mm, ~245 pcs/oz. 1.5-18	\$250.00	—
22 mm, ~164 pcs/oz. 1.5-22	—	\$223.00
25 mm, ~127 pcs/oz. 1.5-25	—	\$252.00
Circles, Erie, USA 	1 oz	cs/10 oz
40 mm, ~54 pcs/oz. 1.5-40	\$29.00	244.00
40 mm required for Biotech's FCS2 chamber.		

*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Counting Chamber - see Counting Chamber Coverslips in
MICROSCOPE Accessories & Parts section

Fused Quartz, **USA**

0.15-0.25 mm Thick, NSG Type OZ®	
1" Square, each. QCS1X1	\$19.00
1" Circle, each. QCS1X0	\$19.00
0.16 ± 0.01 mm Thick, A1	
25 mm Square, each. R525025	\$32.40
25 mm Circle, each. R525000	\$23.40

Fused Silica, **USA**

0.15-0.25 mm Thick, Corning 7980	
1" Square, each. FSCS1X1	\$21.50
1" Circle, each. FSCS1X0	\$21.50
0.16 ± 0.01 mm Thick, S1-UV	
25 mm Square, each. R425025	\$32.00
25 mm Circle, each. R425000	\$28.80