

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker. Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot® , permanent, extra fine point	<u>ea</u>	<u>bx/12</u>
Black. 44102	\$1.70	17.45
Blue. 44103	\$1.70	17.45
Green. 44105	\$1.70	17.45
Red. 44104	—	17.45



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Pens, Super^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic		
Regular Tip (4 mm tapered nib). K-600	\$59.00	
(71310, 195505)		
Fine Tip (2.5 mm tapered nib). K-650	\$44.00	
(71312, 195506)		

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA		
China Marker (Wax Pencil), Bx/12. 02089	\$11.50	

Ink resists xylene, alcohol, acetone and formalin. Remove with acetic acid, chloroform, or ethanol/xylene mixture.



Securiline® Marker II® , permanent, PDC	<u>ea</u>	<u>bx/10</u>
Black, 0.6 mm tip. 1450-20-PDC	\$2.50	20.30
(14-905-30, 72107)		



Sharpie Retractable , permanent		
Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 36701	—	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 1735790	\$2.85	29.65
Blue. 1735792	—	29.65

Continued




Scribes , retractable		
Carbide-Tipped (resharpenable). RSC	\$15.30	
Carbide Tip (replacement). RCT-1	\$6.20	
Diamond-Tipped. RSD	\$22.75	
Diamond Tip (replacement). RDT-1	\$12.90	



Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers		
Small, 26 x 18 cm (10 x 7"). XH-2002	\$251.00	
Large, 63 x 20 cm (25 x 8"). XH-2001	\$335.00	

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers		
10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00



20-Slide, Simport		
Clear Lid (PETG), black base (ABS), each. M920-1	\$160.00	
Black Lid (ABS), black base (ABS), each. M920-2	\$160.00	

Continued