




Staining Dishes

Glass (flint)

- 20-Slide, holds slides horizontally in metal rack
 Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)
 Set (dish, rack & lid). **SDS-20** \$26.60
 Dish & Lid. **SD-20** \$16.20
 Rack only, stainless steel. **SR-20** \$13.30
- 30-Slide, holds slides vertically in metal rack
 Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)
 Set (dish, rack & lid). **SDS-30** \$36.35
 Dish & Lid. **SD-30** \$18.05
 Rack only, stainless steel. **SR-30** \$22.25
- 60-Slide, holds slides vertically in metal rack
 Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)
 Set (dish, rack & lid). **SDS-60** \$49.90
 Rack only, stainless steel. **SR-60** \$28.15



Plastic Set

Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

- Plastic (polyoxymethylene [POM]), holds 25 slides
 Set (dish, rack & lid). **0990-99** \$50.90



674



600P90101



675

Staining Jars, hold slides (up to 76 x 26 mm) vertically

- Coplin, holds 10 slides back-to-back ea cs/12
 Glass w/Ground Glass Cover, USA. **673** \$11.00 122.00
 Glass Cover, Fits 673 Jar. **673-LID** \$4.20 —
- ea cs/12
 Glass w/Screw-On Phenolic Cap, USA. **674** \$11.00 122.00
 Phenolic Cap, Fits 674 Jar, black. **674-LID** ...\$4.25 —
- ea cs/12
 Polypropylene w/Screw-On Cap. **600P90101** \$4.60 45.35
 Plastic jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is flat.
- ea cs/12
 Wide Top, Glass w/Square Lid, USA. **675** \$14.80 164.00
 holds 16 slides back-to-back
 Glass Lid, Fits 675 Jar. **675-LID** \$4.55 —




Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Staining Jars (cont)

Wide Top, Plastic (POM)

- 12-Slide, EasyDip™, Simport ea cs/6
 Jar only, Blue. **M900-12B** \$10.50 57.60
 Jar only, Green. **M900-12G** \$10.50 57.60
 Jar only, Pink. **M900-12P** \$10.50 57.60
 Jar only, White. **M900-12W** \$10.50 57.60
 Jar only, Yellow. **M900-12Y** \$10.50 57.60
 Rack/Lid only, gray. **M905-12DGY** \$11.65 64.00



Anodized aluminum staining jar holder

(Jars sold separately)

- Staining Jar Holder, 6-Place. **M906** \$66.65



4501030



4501042

- 25-Slide, light green, microwaveable, acetal
 Jar & Lid. **4501030** \$11.35
 Rack only. **4501042** \$18.90
 Staining Jar Station, 12-Place. **4501011** \$259.00
 Includes metal frame, (12) 25-slide dishes w/ lids (1030) & (1) rack (1042).
 Dimensions: 29(L) x 5(W) x 5"(H)