



Tissue (Embedding) & Biopsy Cassettes



Base Molds, plastic, disposable, Simport

- 15 x 15 x 5 mm, Cs/1000. **M475-2** \$128.00
- 24 x 24 x 5 mm, Cs/1000. **M475-3** \$128.00
- ★ 30 x 24 x 5 mm, Cs/1000. **M475-4** \$128.00
- 37 x 24 x 5 mm, Cs/1000. **M475-5** \$128.00
- 37 x 24 x 10 mm, Cs/500. **M475-10** \$66.50

Designed for Simport Macrosette™.




Biopsy Bags

- Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011** \$235.00
Solvent resistant, 0.2 x 0.2 mm mesh.
- Paper (Tea Bags), 2¼ x 2¼", Pk/500. **56100** \$38.00

M476-1


Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



Biopsy Pads

- Rectangular
 - Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014** \$40.00
 - Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1** \$188.00
- Circular, Simport
 - Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2** \$228.00


Biopsy Cassettes




Economy Biopsy Cassettes

Biopsy Cassettes, Detached Lids (packed separately)

Economy, White, small square pores, Pk/500. **BC-0109W** \$42.90



Simport Histosette I



Simport Micromesh

Biopsy Cassettes, Attached Lids, acetyl polymer

Histosette™ I, 1 mm square pores, Pk/1500 1 pk pk/3
\$124.00 351.00
30° front labeling area.


- Blue. **M491-6** Pink. **M491-3**
- Green. **M491-4** White. **M491-2**

Uniset™, hinged lids, Pk/1500 1 pk pk/3
\$132.00 379.00
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

- Blue. **M506-6** Pink. **M506-3**
- Green. **M506-4** White. **M506-2**

Micromesh™, 0.38 mm pores, Pk/1000 1 pk pk/3
\$209.00 583.00
45° front labeling area.

- Blue. **M507-6** Pink. **M507-3** Yellow. **M507-5**
- Green. **M507-4** White. **M507-2**



Compatible with Primera Signature™ and Tissue-Tek SmartWrite™ cassette printers.

Biopsy Cassettes, Taped


Micromesh™, 0.38 mm pores, Pk/2000 \$349.00
45° front labeling area.

- Blue. **M407-6T** Pink. **M407-3T** Yellow. **M407-5T**
- Green. **M407-4T** White. **M407-2T**

Continued

Tissue Cassettes

Economy
Tissue
Cassettes



Tissue Cassettes, Detached Lids, slotted pores

Economy, Pk/500 \$42.90
 Blue. **EC-0109B** Pink. **EC-0109P**
 Green. **EC-0109G** White. **EC-0109W**
 Lilac. **EC-0109L**

Histosette™ II, Simport, Pk/1000 \$85.00 1 pk pk/3
 45° front labeling area. 245.00

Blue. **M485-6** Pink. **M485-3**
 Green. **M485-4** White. **M485-2**

Uni-Cassettes™, **USA**, Pk/500 \$59.85

Aqua. **6059** Lime. **6062**
 Blue. **6054** Peach. **6064**
 Burnt Orange. **6056** Pink. **6052**
 Chartreuse. **6061** Red. **6060**
 Gold. **6058** Rose. **6063**
 Green. **6053** Tan. **6055**
 Lavender. **6065** White. **6050**
 Lilac. **6057** Yellow. **6051**



Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores

Histosette™ I, snap-apart lids, Pk/1500 \$124.00 1 pk pk/3
 30° front labeling area. 351.00

Blue. **M490-6** Pink. **M490-3**
 Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500 \$166.00 1 pk pk/3
 45° front labeling area. 477.00

Blue. **M492-6** Pink. **M492-3**
 Green. **M492-4** White. **M492-2**

Uniset™, hinged lids, Pk/1500 \$133.00 1 pk pk/3
 35° front labeling area. (Alternative to Sakura's Uni-Cassette®) 379.00

Blue. **M505-6** Pink. **M505-3**
 Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** \$175.00
 Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all
available metal lids.

Tissue Cassettes, without Lids, acetyl polymer, premium

Capsettes™, **USA**, Pk/500 \$41.20

Aqua. **6039** Lime. **6042**
 Blue. **6034** Peach. **6044**
 Burnt Orange. **6036** Pink. **6032**
 Chartreuse. **6041** Red. **6040**
 Gold. **6038** Rose. **6043**
 Green. **6033** Tan. **6035**
 Lavender. **6045** White. **6030**
 Lilac. **6037** Yellow. **6031**

Simport, 45° front labeling area

Green, Pk/1500. **M480-4** \$93.30

Cassette Storage

Box, Cardboard, white

14 x 14 x 2", Pk/50. **8342** \$121.00

Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
in 4 drawers, each divided into 12 rows.
Cabinets can be stacked up to 8 high.
Dimensions: 49.5 x 38 x 21 cm.



Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$239.00



Simport

Store 1500 cassettes or 990 embedding
rings in 6 stackable drawers.
High-impact polystyrene resists
temperatures from -80°C to +80°C.
Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6 \$152.00

Ultra pure and nonflammable;
chills to -62°F.

Freeze Spray, 10 oz Can. 49185 \$12.00



Continued