


**Tissue Cassettes**

Economy  
Tissue  
Cassettes



**Tissue Cassettes, Detached Lids, slotted pores**

Economy, Pk/500 ..... \$42.90  
 Blue. **EC-0109B**                      Pink. **EC-0109P**  
 Green. **EC-0109G**                      White. **EC-0109W**  
 Lilac. **EC-0109L**

Histosette™ II, Simport, Pk/1000 ..... \$85.00      1 pk      pk/3  
 45° front labeling area.                      245.00  
 Blue. **M485-6**                      Pink. **M485-3**  
 Green. **M485-4**                      White. **M485-2**

Uni-Cassettes™, **USA**, Pk/500 ..... \$59.85  
 Aqua. **6059**                      Lime. **6062**  
 Blue. **6054**                      Peach. **6064**  
 Burnt Orange. **6056**                      Pink. **6052**  
 Chartreuse. **6061**                      Red. **6060**  
 Gold. **6058**                      Rose. **6063**  
 Green. **6053**                      Tan. **6055**  
 Lavender. **6065**                      White. **6050**  
 Lilac. **6057**                      Yellow. **6051**



Simport



Simport  
Macrosette

**Tissue Cassettes, Attached Lids, slotted pores**

Histosette™ I, snap-apart lids, Pk/1500 ..... \$124.00      1 pk      pk/3  
 30° front labeling area.                      351.00  
 Blue. **M490-6**                      Pink. **M490-3**  
 Green. **M490-4**                      White. **M490-2**

Histosette™ II, integral lids, Pk/1500 ..... \$166.00      1 pk      pk/3  
 45° front labeling area.                      477.00  
 Blue. **M492-6**                      Pink. **M492-3**  
 Green. **M492-4**                      White. **M492-2**

Uniset™, hinged lids, Pk/1500 ..... \$133.00      1 pk      pk/3  
 35° front labeling area. (Alternative to Sakura's Uni-Cassette®)      379.00  
 Blue. **M505-6**                      Pink. **M505-3**  
 Green. **M505-4**                      White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** ..... \$175.00  
 Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all  
available metal lids.

**Tissue Cassettes, without Lids, acetyl polymer, premium**

Capsettes™, **USA**, Pk/500 ..... \$41.20  
 Aqua. **6039**                      Lime. **6042**  
 Blue. **6034**                      Peach. **6044**  
 Burnt Orange. **6036**                      Pink. **6032**  
 Chartreuse. **6041**                      Red. **6040**  
 Gold. **6038**                      Rose. **6043**  
 Green. **6033**                      Tan. **6035**  
 Lavender. **6045**                      White. **6030**  
 Lilac. **6037**                      Yellow. **6031**

Simport, 45° front labeling area

Green, Pk/1500. **M480-4** ..... \$93.30

**Cassette Storage**

**Box, Cardboard, white**

14 x 14 x 2", Pk/50. **8342** ..... \$121.00  
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
in 4 drawers, each divided into 12 rows.  
Cabinets can be stacked up to 8 high.  
Dimensions: 49.5 x 38 x 21 cm.



**Tissue Cassette Cabinet, Steel, 4-drawer. 485400** ..... \$239.00



Simport

Store 1500 cassettes or 990 embedding  
rings in 6 stackable drawers.  
High-impact polystyrene resists  
temperatures from -80°C to +80°C.  
Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6** ..... \$152.00

Ultra pure and nonflammable;  
chills to -62°F.

**Freeze Spray, 10 oz Can. 49185** ..... \$12.00



Continued