

# GENERAL LAB SUPPLIES & EQUIPMENT


**Alcohol Prep Pads, Sterile, 70% Isopropanol**  
 Large, 4.5 x 9 cm (1.8 x 3.5"), Cs/1000. **44058-404**..... \$40.00



**Applicators, 6" long**  
 Cotton Tipped, wooden shaft  
 Nonsterile  
 Bg/100. **95-8702** ..... \$1.50  
 Bx/1000, 10 bg/bx. **95-8702** ..... \$12.00  
 Sterile Pk/1, Bx/100. **4404304** ..... \$7.80  
 Sterile Pk/2, Bx/200. **44024-106-2S** ..... \$7.90  
 Plain, wood, Bg/1000. **95-8701** ..... \$7.25



**Aprons**  
 Disposable, polyethylene, 46"(L) x 28"(W)  
 0.8 mil Thick (0.02 mm), Bx/100. **105DA2846** ..... \$10.95  
 1.75 mil Thick (0.044 mm), Bx/100. **59010485** ..... \$38.00  
 Rubberized Cloth, black  
 36"(L) x 27"(W), Pk/10. **3101276010** ..... \$66.50  
 42"(L) x 27"(W), Pk/10. **3101276110** ..... \$76.00  
 ★ 42"(L) x 36"(W), Pk/10. **3101276310** ..... \$90.00  
 46"(L) x 36"(W), Pk/10. **3101276210** ..... \$96.30




**Autoclaves, benchtop, Benchmark (1-year warranty)**  
 8 L, Bioclave™ Mini. **B4000-M** ..... \$3,219.00  
 16 L, Bioclave™. **B4000-16** ..... \$4,377.00

## Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



**Bag, Biohazard, Benchtop, Bx/500. HS10320** ..... \$69.35  
**Biohazard Bag Holder, Benchtop. HS10044H**..... \$20.85  
 Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

**Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick**  
 14 x 19" (36 x 48 cm), Bx/200. **131641419** ..... \$69.00  
 19 x 23" (48 x 58 cm), Bx/200. **131641923** ..... \$99.00  
 25 x 35" (64 x 89 cm), Bx/200. **131642535** ..... \$229.00

**Bags, Biohazard, Heavy-Duty, 2 mil (0.05 mm) thick**  
 14 x 19" (36 x 48 cm), Bx/200. **131651419** ..... \$191.00  
 19 x 23" (48 x 58 cm), Bx/200. **131651923** ..... \$279.00  
 25 x 35" (64 x 89 cm), Bx/200. **131652535** ..... \$505.00



**Bags, Reclosable, polyethylene**  
 2 mil Thick, Clear  
 3 x 7", Pk/100. **17730-100** ..... \$7.20  
 2 mil Thick, White Block for labeling  
 3 x 4", Pk/100. **2012-100** ..... \$7.40  
 6 x 9", Pk/100. **1309-100** ..... \$9.20  
 6 x 12", Pk/100. **10850-100** ..... \$15.40  
 4 mil Thick, clear  
 3 x 4", Pk/100. **494-100** ..... \$12.00

Safely transport patient specimens with these OSHA and NCCLS compliant bags. Each resealable, 2 mil thick polyethylene bag has an attached document pouch and imprinted biohazard symbol.



**Bag, Specimen Transport, 6 x 9", Bx/1000. 4919**..... \$56.70



**Bags, Waste Disposal, HDPE, Cleanware™**  
 8 x 10", Pk/50. **131741008** ..... \$19.90  
 12 x 16", Pk/50. **131741012** ..... \$33.45