



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

- Bag Sealer**, electric, 12" (30.5 cm) Wide. **HS28272**..... \$168.00
Bag Sealer Replacement Part Kit. HS28280..... \$30.65
Poly Tubing
 4 mil
 4" x 1500' Roll. **1143** \$110.00
 6" x 1500' Roll. **1144** \$162.00
 6 mil
 4" x 1000' Roll. **1658** \$132.00
 6" x 1000' Roll. **1659** \$162.00
Poly Tubing Dispenser. H-28..... \$34.00

Balances

Balance, Portable



Model	DWP-1001
Capacity (g)	1000
Readability (g)	0.1
Modes	g, oz, ct, gn
Pan	14.6 cm (5.75") stainless steel
Power	4 AA/AC adapter
Price	\$72.00

Warranty: 1 year

Benchmark
W3200-1200



Balances, Toploading

Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$999.00	\$808.00

Warranty: 2 years

Benchmark



Balances, Analytical

Model	W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,487.00	\$1,580.00

Warranty: 2 years

Baths



Dry Baths, digital, Benchmark (2-year warranty)

- 1-Block Capacity. **BSH1001**..... \$359.00
 2-Block Capacity. **BSH1002**..... \$425.00
 4-Block Capacity. **BSH1004**..... \$725.00

Blocks sold separately – please specify configuration

Blocks, aluminum

- Holds (48) 0.2 mL Tubes. **BSW02**..... \$85.00
 Holds (24) 0.5 mL Tubes. **BSW05**..... \$85.00
 Holds (24) 1.5 mL Tubes. **BSW1500**..... \$85.00
 Holds (24) 1.5/2.0 mL Tubes. **BSW1520**..... \$85.00
 Holds (20) 12/13 mL Tubes. **BSW13**..... \$85.00
 Holds (12) 15 mL Tubes. **BSW15**..... \$85.00
 Holds (5) 50 mL Tubes. **BSW50**..... \$85.00



Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

- 4 L, 3 tube racks. **B2000-4** \$739.00
 8 L. **B2000-8**..... \$799.00
 12 L. **B2000-12**..... \$956.00

Continued