

## Carboys & Jugs



### Carboys with Spigot, polypropylene, round

10 L (2.5 gal). <b>118460025</b> .....	\$55.50
20 L (5 gal). <b>60034052</b> .....	\$72.30



### Jug, dispensing, HDPE jug/LDPE spigot

20 L (5 gal). <b>118500000</b> .....	\$46.75
--------------------------------------	---------



### Carts, utility, stainless steel, Lakeside

200 lb Capacity

27½"(L) x 16¼"(W) x 32⅞"(H). <b>311</b> .....	\$272.00
---	----------

*Model 311 can ship UPS*

30¾"(L) x 18⅝"(W) x 33"(H). <b>322</b> .....	\$311.00
--	----------



C1012



C1008-C



MLX-110

### Centrifuge, Mini

5,500 RPM (2000 x g), 12-place, MyFuge™ 12. <b>C1012</b> .....	\$321.00
--	----------

1-year warranty

6,000 RPM (2000 x g), 8-place, MyFuge™ Mini. <b>C1008-C</b> .....	\$229.00
---	----------

1-year warranty

10,000 RPM (5000 x g), 6-place, blue body. <b>MLX-110</b> .....	\$185.00
---	----------

2-year warranty



## Clamps & Supports

Rectangular Steel Base



Triangular Cast-Iron Base



### Support (Ring) Stands, finished metal base/plated support rod

#### Rectangular Steel Base

6 x 4" base, 5/16 x 18" rod. <b>3104SSS66</b> .....	\$12.80
---	---------

8 x 5" base, 3/8 x 20" rod. <b>3104SSS67</b> .....	\$15.85
--	---------

9 x 6" base, 1/2 x 24" rod. <b>3104SSS68</b> .....	\$18.55
--	---------

#### Rectangular Cast-Iron Base

Provides greater stability than steel base.

6 x 4" base, 5/16 x 18" rod. <b>3104SSC14</b> .....	\$14.35
---	---------

8 x 5" base, 3/8 x 20" rod. <b>3104SSC15</b> .....	\$16.90
--	---------

9 x 6" base, 1/2 x 24" rod. <b>3104SSC16</b> .....	\$23.00
--	---------

#### Triangular Cast-Iron Base

5" Legs, 1/2 x 24" rod. <b>1257G2324</b> .....	\$26.85
--	---------

6" Legs, 1/2 x 36" rod. <b>1257G23</b> .....	\$30.65
--	---------



3104BC37




600THCL01



600CLBURDD



97-4023



1257G49

### Support Clamps

Burette, plastic coated round jaws, 360° rotatable. <b>3104BC37</b> .....	\$6.90
---	--------

Burette, Double, deluxe. <b>600CLBURDD</b> .....	\$18.30
--	---------

Plastic covered spring jaws on a metal frame securely grip burettes up to 30 mm OD. Fits rods up to ¾" (19 mm) OD.

Hook Collar, nickel plated. <b>1257G49</b> .....	\$5.45
--	--------

Fits rods up to 1/2" (13 mm) in diameter.

Right-Angle, aluminum. <b>97-4023</b> .....	\$6.10
---	--------

Fits rods up to 5/8" in diameter.

Thermometer, aluminum, rotates 360°. <b>600THCL01</b> .....	\$7.45
---	--------

Fits rods up to 1/2" (13 mm) in diameter.



### Support Rings, plated cast-iron

3" OD (2½" ID). <b>600SRC103</b> .....	\$5.70
--	--------

4" OD (3½" ID). <b>600SRC104</b> .....	\$7.40
--	--------

5" OD (4½" ID). <b>600SRC105</b> .....	\$7.85
--	--------