



Non-Vacuum

Glazed porcelain sample plate has 3.5-4.5 cm center hole.

Vacuum



Desiccators

- Non-Vacuum, glass, economy
24 cm (10") ID. **125720002**.....\$65.75
- Vacuum, polycarbonate top/polypropylene bottom
138 mm Plate Diameter. **60055204**.....\$52.45
- 190 mm Plate Diameter. **60055205**.....\$71.25
- 238 mm Plate Diameter. **60055206**.....\$90.30



- Dippers**, 500 mL polypropylene cup with HDPE handle
- 3' Handle. **367840016**\$52.15
 - 6' Handle. **367850016**\$67.95
 - Replacement Cup, screw-on. **367860016**\$12.40



Autoclavable polypropylene bottle with square bottom holds 170 mL (6 oz). Simport

Drosophila Stock Bottle, Cs/500. **B352**.....\$142.00

Drosophila Vials *Call us to inquire or get a quote*



188050000



HS23243A



188190023

Drying Racks

- Non-Electric
Bel-Art. **188050000**.....\$103.00
Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 17 x 43 x 48 cm(H).
- Heathrow. **HS23243A**.....\$110.00
- Electric, double-sided, 3 tier, Bel-Art. **188190023**.....\$1,025.00
1-year warranty
- Dust Filter & Holder (replacement). **999321563**.....\$12.00
- Heating Element (replacement). **999321564**.....\$193.00



Eyewash Station, single 32 oz bottle (empty). **248660000**.....\$47.30

Eyewear

Absorbs 99.9% UV



Glasses, Safety, integral side shield. **43574002**.....\$5.60



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**.....\$3.70

Film & Foil



- Parafilm**, thermoplastic sealing film, **USA**
- 2" x 250' (5 cm x 76 m) Roll. **HS234526A**.....\$26.50
 - 4" x 125' (10 cm x 38 m) Roll. **HS234526B**.....\$27.00
 - 4" x 250' (10 cm x 76 m) Roll. **HS234526C**.....\$45.20

Holds one or two 2" rolls, or one 4" roll.



Parafilm Dispenser, ABS plastic, blue. **HS234525B**.....\$38.70



- Foil**, Super-Polyfoil™, polyester-backed aluminum foil
- 4 x 4" (10.2 x 10.2 cm) Sheets, Pk/100. **132000044**.....\$34.65
 - 5 1/8 x 5 1/8" (13 x 13 cm) Sheets, Pk/100. **132000066**.....\$44.95