


# pH Meters, Hanna



Model	HI98103	HI98128	HI8424	HI2221	HI2210
<b>Modes</b>	pH	pH, °C/°F	pH, mV (*ORP), °C/°F	pH, ORP, °C/°F	pH, °C/°F
<b>pH</b>	<b>Range</b>	0.00 to 14.00 pH	-2.00 to 16.00 pH	-2.00 to 16.00 pH	-2.00 to 16.00 pH
	<b>Resolution</b>	0.01 pH	0.01 pH	0.01 pH	0.01 pH
	<b>Accuracy</b>	±0.2 pH	±0.05 pH	±0.01 pH	±0.01 pH
<b>mV</b>	<b>Range</b>	-	-	±699.9; ±1999 mV	±699.9; ±2000 mV
	<b>Resolution</b>	-	-	0.1; 1 mV	0.1; 1 mV
	<b>Accuracy</b>	-	-	±0.2; ±1 mV	±0.2; ±1 mV
<b>Temp.</b>	<b>Range</b>	-	-5.0 to 60.0°C	-20.0 to 120.0°C	-20.0 to 120.0°C
	<b>Resolution</b>	-	0.1°C	0.1°C	0.1°C
	<b>Accuracy</b>	-	±0.5°C	±0.4°C	±0.2°C
<b>Calibration Points</b>	2	2	2	A/M	5
<b>Temp. Compensation: Auto (A) or Manual (M)</b>	-	A	A	12V DC adapter	A/M
<b>Power</b>	AAA (2) batteries	AAA (4) batteries	9V battery		12V DC adapter
<b>Features</b>	<ul style="list-style-type: none"> <li>• Pocket-size</li> <li>• Replaceable electrode</li> </ul>	<ul style="list-style-type: none"> <li>• Waterproof</li> <li>• Replaceable electrode</li> </ul>	<ul style="list-style-type: none"> <li>• Double-junction, gel-filled electrode</li> </ul>	<ul style="list-style-type: none"> <li>• Calibration Check™ technology</li> <li>• On-screen electrode condition and response time</li> <li>• GLP features</li> <li>• USB port</li> <li>• Electrode holder</li> <li>• Refillable electrode</li> </ul>	<ul style="list-style-type: none"> <li>• Electrode holder</li> <li>• Refillable electrode</li> </ul>
<b>Price</b>	\$40.60	\$106.00	\$333.00	\$488.00	\$425.00

\*Additional ORP (Oxidation Reduction Potential) electrode required


Warranty: 2 years on meter; 1 year on electrode



Measure pH with 1 pH unit accuracy from 1-14. Package contains 100 paper strips, 2 x 1/4", in a plastic vial with color chart.

**pH Paper, pH 1-14, Pk/100. 600PHP100 ..... \$3.70**

### Pipetting



- Dropper Bulbs**
- Amber Latex
    - 1 mL Capacity, 43 mm long, Pk/12. **20621A** ..... \$6.20
    - 1 mL Capacity, 43 mm long, Pk/72. **HS20622A** ..... \$32.55
    - 2 mL Capacity, 52 mm long, Pk/12. **20621B** ..... \$6.20
    - 2 mL Capacity, 52 mm long, Pk/72. **HS20622B** ..... \$32.55
  - Natural Rubber
    - 1 mL Capacity, 30 mm long, black, Pk/12. **200-1B** ..... \$3.50
    - 2 mL Capacity, 32 mm long, red, Pk/12. **200-2B** ..... \$5.85
  - Natural Rubber w/Flange, blue
    - 1 mL Capacity, 29 mm long, Pk/24. **200-1C** ..... \$12.15

*Continued*