


Medicine (Eye) Dropper, Pk/12. 97-4301 \$3.60



Pipet Fillers

3-Valve, rubber, black. **207B** \$14.40
 Wheel-Type, Pipette Pump™
 Fits Pipets Up To 2 mL, blue. **378970000** \$24.35
 Fits Pipets Up To 10 mL, green. **378980000** \$25.45
 Fits Pipets Up To 25 mL, red. **378990000** \$35.05
 Wheel-Type, Pipette Pump™ II, quick-release ea pk/12
 Fits Pipets Up To 2 mL, blue. **379111002** \$26.35 299.00
 Fits Pipets Up To 10 mL, green. **379111010** \$26.35 299.00
 Fits Pipets Up To 25 mL, red. **379111025** \$35.05 397.00



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, blue, Heathrow. HS20612M \$72.45

Pipets, Mohr (Measuring), Class A, individually certified


5 mL, (±0.020 mL), Pk/5. **600PM7058H** \$64.80
 10 mL, (±0.030 mL), Pk/5. **600PM7058I** \$64.80
 25 mL, (±0.050 mL), Pk/5. **600PM7058J** \$96.60



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™

Nonsterile
 5 mL, blue band, pk/25, Cs/375. **1755** \$74.40
 10 mL, orange band, pk/25, Cs/250. **1775** \$50.85
 25 mL, red band, pk/10, Cs/100. **1784** \$47.35
 Sterile Pk/1 (individually paper/plastic wrapped)
 1 mL, yellow band, Cs/500. **1700** \$74.65
 2 mL, green band, Cs/500. **1720** \$80.05
 5 mL, blue band, reverse grad., Cs/250. **1740** \$59.05
 10 mL, orange band, reverse grad., Cs/200. **1760** \$51.50
 25 mL, red band, reverse grad., Cs/100. **1780** \$55.65
 50 mL, black band, reverse grad., Cs/100. **1820** \$139.00



Pipets, Transfer, LDPE (low-density polyethylene)

Graduated, Nonsterile, 15 cm long
 5 mL capacity, 1 mL graduated, 25 drops/mL
 Bx/500. **PTP-02** \$16.45
 Cs/10 Bx. **PTP-02** \$153.00
 7 mL capacity, 3 mL graduated, 14 drops/mL
 Bx/500. **PTP-01** \$16.45
 Cs/10 Bx. **PTP-01** \$153.00
 Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
 5 mL capacity, 1 mL graduated, Simport
 Bx/400. **P200-521S** \$39.65
 Cs/10 Bx. **P200-521S** \$364.00
 Graduated, Sterile Pk/20, 15 cm long
 5 mL capacity, 1 mL graduated, Bx/500. **HS206371B** \$33.30
 Long Stem, Sterile Pk/1 (individually wrapped)
 22.5 cm (9" long), 6 mL capacity, Cs/400. **139038** \$59.45
 Long Stem, Nonsterile
 30 cm (12" long), 23 mL capacity, Bx/100. **139050** \$17.70

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



Pipets, Volumetric, Class A, individually certified

2 mL (±0.006 mL), orange band, Pk/5. **600PT7100B** \$65.30
 5 mL (±0.010 mL), white band, Pk/5. **600PT7100C** \$68.65
 10 mL (±0.020 mL), red band, Pk/5. **600PT7100D** \$78.40
 20 mL (±0.030 mL), yellow band, Pk/5. **600PT7100E** \$111.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond Jr.
 5 µL, red. **3322-5** \$18.20
 10 µL, pink. **3322-10** \$18.20
 20 µL, green. **3322-20** \$18.20
 50 µL, purple. **3322-50** \$18.20
 100 µL, orange. **3322-100** \$18.20



Adjustable Volume, Diamond Advance™ (3-year warranty)
 0.5-10 µL, 0.02 µL increments, red. **3341-10** \$219.00
 5-50 µL, 0.1 µL increments, yellow. **3341-50** \$219.00
 10-100 µL, 0.1 µL increments, blue. **3341-100** \$219.00
 100-1000 µL, 1 µL increments, orange. **3341-1000** \$219.00

Continued