


Pipettor Holders

Carousel, Flip & Grip™, Heathrow

- Blue. **HS120026** \$76.20
 Lime. **HS120056** \$76.20



Stands, Heathrow


- 4-Place, clear acrylic, **USA. HS206204** \$39.30
 6-Place
 Blue, opaque ABS plastic. **HS20613E** \$38.25
 Clear, acrylic, **USA. HS20620** \$50.90
 (82024-544)



Pipettor Tips, Filtered (Barrier) Access www.globescientific.com
 and call us to place an order

Pipettor Tips, Standard, Sterile, graduated

- 100-1250 µL, blue, 6 racks of 96, Bx/576. **151153BRS-96** \$38.55



Reagent Reservoirs

- 25 mL, white polystyrene
 Sterile Bg/5, Bx/200. **P8025-5S** \$71.30
 55 mL, v-bottom, Heathrow
 Nonsterile, clear PVC, Pk/100. **HS20521A** \$40.20
 Sterile Bg/1, white polystyrene, Pk/50. **HS20521C** \$69.80



Polypropylene




Stainless Steel

Pitchers

- Polypropylene, graduated (mL), tall form, tapered
 1000 mL, 50 mL graduations, each. **289910000** \$6.10
 2000 mL, 100 mL graduations, each. **289920000** \$6.90
 3000 mL, 100 mL graduations, each. **289930000** \$13.80
 Stainless Steel, graduated (mL/oz), tapered, Polar Ware
 1000 mL/32 oz, each. **T1063** \$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

- 2-D (seesaw), 12 x 12" platform w/flat mat
 BenchBlotter™. **BR1000** \$413.00
 • Fixed speed of 12 rocks/minute
 • Platform angle is ±8°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 2.3 kg (5 lbs)
 3-D (nutating)
 Small, 10.5 x 7.5" platform w/dimpled mat
 MiniMixer™. **B3D1020** \$350.00
 Large, 12 x 12" platform w/dimpled mat
 BioMixer™. **B3D1320** \$422.00
 • Fixed speed of 24 rocks/minute
 • Platform angle is ±20°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 1.8 kg (4 lbs)



BR2000



B3D2300

Rockers, Variable Speed

- Benchmark (2-year warranty)
 2-D (seesaw), 14 x 12" platform w/flat mat
 BenchRocker™ 2D. **BR2000** \$679.00
 Stacking Platform. **BR2000-STACK** \$102.00
 14 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP** \$36.10
 Adds 1.25" to platform separation.
 • Variable speed of 2-30 rocks/minute
 • Platform tilt angle adjusts from 0-30°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 3.6 kg (8 lbs)
 3-D (nutating), 12 x 12" platform w/dimpled and flat mats
 BenchRocker™ 3D. **B3D2300** \$722.00
 Stacking Platform. **BR1000-STACK** \$94.00
 12 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP** \$36.10
 Adds 1.25" to platform separation.
 • Variable speed of 2-30 rocks/minute
 • Platform tilt angle adjusts from 0-30°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 2.3 kg (5 lbs)

