


## Rulers, Meter Sticks & Tape Measure

Stainless Steel



**Rulers**, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461** ..... \$0.95

White. **4460** ..... \$0.55

Stainless Steel, 12 mm wide, **USA. 605-065** ..... \$5.50

**Meter Sticks**, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. **605-050** ..... \$4.90

Metal Ends. **605-051** ..... \$6.15

**Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036** ..... \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions.

The opposite side is graduated in inches with 1/16" divisions.



440065



4402269

**Sharps Container**, red polypropylene

1 Quart, translucent locking lid. **440065** ..... \$3.80

Dimensions: 4.25" x 4.25" x 6.25"(H)

5.4 L, horizontal entry lid, Prevent™, each. **4402269** ..... \$7.95

Dimensions: 12" x 4.75" x 11"(H)

Screw-Cap, 90 mL



**Specimen Containers**, polypropylene (PP)

Screw-Cap

90 mL, sterile, Cs/400. **840175** ..... \$110.00

480 mL (16 oz), 83 mm (3 1/4") opening, Cs/24. **14510** ..... \$71.95

960 mL (32 oz), 111 mm (4 3/8") opening, Cs/12. **14511** ..... \$59.05

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. **271016W** ..... \$57.15

9.9 cm opening x 11.1 cm high

1000 mL (34 oz), Cs/100. **271032W** ..... \$83.75

11.7 cm opening x 13.7 cm high.



10"(D) x 12"(W) x 18"(H)

**Splash Shield**, 3/16" acrylic, angled top, TrippNT. **50707** ..... \$49.20



**Spot Plates**, glazed white porcelain, 21 mm diameter wells

6-Well. **LPC-161** ..... \$2.35

12-Well, 115 x 95 x 7 mm. **3104187701** ..... \$3.75 ↓

## Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

**Stirring Bars**, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

40 mm (1 1/2") Long. **371160002** ..... \$25.35

50 mm (2") Long. **37116003** ..... \$37.90

Micro, Spinbar®

White

5 x 2 mm. **371190005** ..... \$5.05

7 x 2 mm. **371190007** ..... \$5.05

10 x 3 mm. **371190010** ..... \$5.55

12.7 x 3 mm. **371190127** ..... \$5.55

15 x 1.5 mm. **371190015** ..... \$4.95

Red

7 x 2 mm. **371210013** ..... \$5.55

8 x 1.5 mm. **371210001** ..... \$5.55

Octagonal, molded-on pivot ring, Spinbar®

White

12.7 x 8 mm (1/2 x 5/16"). **371101256** ..... \$4.70

15.9 x 8 mm (5/8 x 5/16"). **371100058** ..... \$4.70

25.4 x 9.5 mm (1 x 3/8"). **371100138** ..... \$6.20

38.1 x 9.5 mm (1 1/2 x 3/8"). **371101128** ..... \$7.40

50.8 x 9.5 mm (2 x 3/8"). **371100238** ..... \$9.90

63.5 x 9.5 mm (2 1/2 x 3/8"). **371102128** ..... \$11.40

76.2 x 12.7 mm (3 x 1/2"). **371100003** ..... \$20.50

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034** ..... \$5.75

Continued