


## Rulers, Meter Sticks & Tape Measure

Stainless Steel



**Rulers**, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461** ..... \$0.95

White. **4460** ..... \$0.55

Stainless Steel, 12 mm wide, **USA. 605-065** ..... \$5.50

**Meter Sticks**, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. **605-050** ..... \$4.90

Metal Ends. **605-051** ..... \$6.15

**Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036** ..... \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



440065



4402269

**Sharps Container**, red polypropylene

1 Quart, translucent locking lid. **440065** ..... \$3.80

Dimensions: 4.25" x 4.25" x 6.25"(H)

5.4 L, horizontal entry lid, Prevent™, each. **4402269** ..... \$7.95

Dimensions: 12" x 4.75" x 11"(H)

Screw-Cap, 90 mL



Screw-Cap,  
480 & 960 mL



Snap-Cap

**Specimen Containers**, polypropylene (PP)

Screw-Cap

90 mL, sterile, Cs/400. **840175** ..... \$110.00

480 mL (16 oz), 83 mm (3 1/4") opening, Cs/24. **14510** ..... \$71.95

960 mL (32 oz), 111 mm (4 3/8") opening, Cs/12. **14511** ..... \$59.05

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. **271016W** ..... \$57.15

9.9 cm opening x 11.1 cm high


1000 mL (34 oz), Cs/100. **271032W** ..... \$83.75

11.7 cm opening x 13.7 cm high.



10"(D) x 12"(W) x 18"(H)

**Splash Shield**, 3/16" acrylic, angled top, TrippNT. **50707** ..... \$49.20



**Spot Plates**, glazed white porcelain, 21 mm diameter wells

6-Well. **LPC-161** ..... \$2.35

12-Well, 115 x 95 x 7 mm. **3104187701** ..... \$3.75 ↓

## Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

**Stirring Bars**, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

40 mm (1 1/2") Long. **371160002** ..... \$25.35

50 mm (2") Long. **37116003** ..... \$37.90

Micro, Spinbar®

White

5 x 2 mm. **371190005** ..... \$5.05

7 x 2 mm. **371190007** ..... \$5.05

10 x 3 mm. **371190010** ..... \$5.55

12.7 x 3 mm. **371190127** ..... \$5.55

15 x 1.5 mm. **371190015** ..... \$4.95

Red

7 x 2 mm. **371210013** ..... \$5.55

8 x 1.5 mm. **371210001** ..... \$5.55

Octagonal, molded-on pivot ring, Spinbar®

White

12.7 x 8 mm (1/2 x 5/16"). **371101256** ..... \$4.70

15.9 x 8 mm (5/8 x 5/16"). **371100058** ..... \$4.70

25.4 x 9.5 mm (1 x 3/8"). **371100138** ..... \$6.20

38.1 x 9.5 mm (1 1/2 x 3/8"). **371101128** ..... \$7.40

50.8 x 9.5 mm (2 x 3/8"). **371100238** ..... \$9.90

63.5 x 9.5 mm (2 1/2 x 3/8"). **371102128** ..... \$11.40

76.2 x 12.7 mm (3 x 1/2"). **371100003** ..... \$20.50

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034** ..... \$5.75

Continued

### Stirring Bars (cont)

- Stirring Bar Assortment, octagonal, Pk/6. **371600000** ..... \$36.55  
Includes two of each size: 25.4 x 8 mm (1 x 5/16"), 38.1 x 8 mm (1 1/2 x 5/16"), and 50.8 x 8 mm (2 x 5/16").
- Polygonal, molded-on pivot ring, Spinbar®  
12 x 4.5 mm (7/16 x 3/16"). **371220012**..... \$5.10
- Triangular, vanes w/cross-mounted magnets, Spinvane®  
Fits 1 mL Vials. **371340000** ..... \$9.85  
Fits 3/5 mL Vials. **371350000** ..... \$9.85



**Stirring Bar Rack, Spinbar® Garage. 377900000** ..... \$46.60



**Stirring Bar Retainer, Spinbar® Restrainer. 377861000**..... \$13.45


### Stirring Bar Retrievers, magnetic, Spinbar®

- 30.5 cm (12") Long. **377720000** ..... \$14.85
- 45.7 cm (18") Long. **377730000** ..... \$17.20



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

**Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000** ..... \$28.95



### Stirring Rods, glass, annealed ends

- 150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** ..... \$4.25
- Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** ..... \$5.95
- 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** ..... \$6.10
- 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012**..... \$22.90



### Stoppers



Stoppers, rubber, black

Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6 1/2	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. <b>50020</b> .....	\$3.10	13.10
Size 3, ~29 pcs/lb. <b>50030</b> .....	\$4.15	13.10
Size 4, ~23 pcs/lb. <b>50040</b> .....	\$5.25	13.10
Size 5, ~20 pcs/lb. <b>50050</b> .....	—	13.10
Size 6, ~15 pcs/lb. <b>50060</b> .....	—	13.10
Size 6.5, ~14 pcs/lb. <b>500650</b> .....	—	13.10
Size 7, ~11 pcs/lb. <b>50070</b> .....	—	13.10
Size 8, ~9 pcs/lb. <b>50080</b> .....	—	13.10
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. <b>50011</b> .....	\$2.35	13.10
Size 2, ~43 pcs/lb. <b>50021</b> .....	\$2.80	13.10
★ Size 3, ~30 pcs/lb. <b>50031</b> .....	\$4.00	13.10
Size 4, ~24 pcs/lb. <b>50041</b> .....	\$5.00	13.10
Size 5, ~21 pcs/lb. <b>50051</b> .....	—	13.10
Size 6, ~18 pcs/lb. <b>50061</b> .....	—	13.10
Size 9, ~8 pcs/lb. <b>50091</b> .....	—	13.10
Size 12, ~3 pcs/lb. <b>50121</b> .....	—	13.10
2-Hole	lb	
Size 6, ~17 pcs/lb. <b>50062</b> .....	\$13.10	
Size 7, ~12 pcs/lb. <b>50072</b> .....	\$13.10	
Size 8, ~9 pcs/lb. <b>50082</b> .....	\$13.10	
Size 9, ~7 pcs/lb. <b>50092</b> .....	\$13.10	
Size 10, ~5 pcs/lb. <b>50102</b> .....	\$13.10	
Assortment, 1 lb. <b>50012</b> .....	\$17.10	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.