

Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461** \$0.95

White. **4460** \$0.55

Stainless Steel, 12 mm wide, **USA. 605-065** \$5.50

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. **605-050** \$4.90

Metal Ends. **605-051** \$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



440065



4402269

Sharps Container, red polypropylene

1 Quart, translucent locking lid. **440065** \$3.80

Dimensions: 4.25" x 4.25" x 6.25"(H)

5.4 L, horizontal entry lid, Prevent™, each. **4402269** \$7.95

Dimensions: 12" x 4.75" x 11"(H)

Screw-Cap, 90 mL



Screw-Cap,
480 & 960 mL



Snap-Cap

Specimen Containers, polypropylene (PP)

Screw-Cap

90 mL, sterile, Cs/400. **840175** \$110.00

480 mL (16 oz), 83 mm (3 1/4") opening, Cs/24. **14510** \$71.95

960 mL (32 oz), 111 mm (4 3/8") opening, Cs/12. **14511** \$59.05

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. **271016W** \$57.15

9.9 cm opening x 11.1 cm high

1000 mL (34 oz), Cs/100. **271032W** \$83.75

11.7 cm opening x 13.7 cm high.



10"(D) x 12"(W) x 18"(H)

Splash Shield, 3/16" acrylic, angled top, TrippNT. **50707** \$49.20



Spot Plates, glazed white porcelain, 21 mm diameter wells

6-Well. **LPC-161** \$2.35

12-Well, 115 x 95 x 7 mm. **3104187701** \$3.75 ↓

Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

40 mm (1 1/2") Long. **371160002** \$25.35

50 mm (2") Long. **37116003** \$37.90

Micro, Spinbar®

White

5 x 2 mm. **371190005** \$5.05

7 x 2 mm. **371190007** \$5.05

10 x 3 mm. **371190010** \$5.55

12.7 x 3 mm. **371190127** \$5.55

15 x 1.5 mm. **371190015** \$4.95

Red

7 x 2 mm. **371210013** \$5.55

8 x 1.5 mm. **371210001** \$5.55

Octagonal, molded-on pivot ring, Spinbar®

White

12.7 x 8 mm (1/2 x 5/16"). **371101256** \$4.70

15.9 x 8 mm (5/8 x 5/16"). **371100058** \$4.70

25.4 x 9.5 mm (1 x 3/8"). **371100138** \$6.20

38.1 x 9.5 mm (1 1/2 x 3/8"). **371101128** \$7.40

50.8 x 9.5 mm (2 x 3/8"). **371100238** \$9.90

63.5 x 9.5 mm (2 1/2 x 3/8"). **371102128** \$11.40

76.2 x 12.7 mm (3 x 1/2"). **371100003** \$20.50

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034** \$5.75

Continued

Novelty Items



Jewelry, sterling silver

Charms

- Caduceus, 7/8" long. **88501** \$11.00
- Chimpanzee, 1 1/16" long. **88503** \$12.00
- Microscope, 3/4" long. **88505** \$11.00
- Mouse, 1" wide. **88507** \$12.00

Earrings, pair

- Stud, microscope, 3/4" long. **88521** \$34.00
- Wire Hoops
 - Microscope, 3/4" long. **88531** \$25.00
 - Mouse, 1" wide. **88533** \$27.00
- Pin, Microscope. **88541** \$17.00

Micro Bible. 89585 \$15.00

View with a 10X objective all 1,245 pages of the King James version of the Holy Bible on this approximately 1" square microfilm.



Ruler, 1 Decimeter, stainless steel, Germany. 605-060 \$10.20

Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.



This beautiful, transparent ornament is a collage of cytological color photomicrographs creatively captured on a 3" glass disk. Suspended from a gold braid.

Suncatcher, "Cytological Crash Course." CCC01 \$22.50



Tube-Wringer®, Heavy-Duty, USA. 295401 \$24.80

3" aluminum rollers with steel handle