

Stirring Bars (cont)

- Stirring Bar Assortment, octagonal, Pk/6. **371600000** \$36.55
Includes two of each size: 25.4 x 8 mm (1 x 5/16"), 38.1 x 8 mm (1 1/2 x 5/16"), and 50.8 x 8 mm (2 x 5/16").
- Polygonal, molded-on pivot ring, Spinbar®
12 x 4.5 mm (7/16 x 3/16"). **371220012**..... \$5.10
- Triangular, vanes w/cross-mounted magnets, Spinvane®
Fits 1 mL Vials. **371340000** \$9.85
Fits 3/5 mL Vials. **371350000** \$9.85



Stirring Bar Rack, Spinbar® Garage. 377900000 \$46.60



Stirring Bar Retainer, Spinbar® Restrainer. 377861000 \$13.45


Stirring Bar Retrievers, magnetic, Spinbar®

- 30.5 cm (12") Long. **377720000** \$14.85
- 45.7 cm (18") Long. **377730000** \$17.20



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000 \$28.95



Stirring Rods, glass, annealed ends

- 150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.25
- Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$5.95
- 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.10
- 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$22.90



Stoppers



Stoppers, rubber, black

Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6 1/2	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.10	13.10
Size 3, ~29 pcs/lb. 50030	\$4.15	13.10
Size 4, ~23 pcs/lb. 50040	\$5.25	13.10
Size 5, ~20 pcs/lb. 50050	—	13.10
Size 6, ~15 pcs/lb. 50060	—	13.10
Size 6.5, ~14 pcs/lb. 500650	—	13.10
Size 7, ~11 pcs/lb. 50070	—	13.10
Size 8, ~9 pcs/lb. 50080	—	13.10
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.35	13.10
Size 2, ~43 pcs/lb. 50021	\$2.80	13.10
★ Size 3, ~30 pcs/lb. 50031	\$4.00	13.10
Size 4, ~24 pcs/lb. 50041	\$5.00	13.10
Size 5, ~21 pcs/lb. 50051	—	13.10
Size 6, ~18 pcs/lb. 50061	—	13.10
Size 9, ~8 pcs/lb. 50091	—	13.10
Size 12, ~3 pcs/lb. 50121	—	13.10
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$13.10	
Size 7, ~12 pcs/lb. 50072	\$13.10	
Size 8, ~9 pcs/lb. 50082	\$13.10	
Size 9, ~7 pcs/lb. 50092	\$13.10	
Size 10, ~5 pcs/lb. 50102	\$13.10	
Assortment, 1 lb. 50012	\$17.10	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.