



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

**Stopper/Tubing Sizer. 133920000** ..... \$13.20

Non-fenestrated towels are sterile, latex-free and poly-lined.



**Surgical Drape, 18 x 26", Bx/50. 44025-515**..... \$25.50

**Tape**



**Tape**

Autoclave, steam


$\frac{3}{4}$ " x 60 Yd Roll. **44073ST048** ..... \$8.50  
 1" x 60 Yd Roll. **44073ST036** ..... \$10.70

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll  
 $\frac{1}{2}$ " Wide, each ..... \$5.35  
 Green. **134623005**    Yellow. **134622005**  
 $\frac{3}{4}$ " Wide, Pk/4 ..... \$23.20  
 Green. **134820075**    Yellow. **134810075**  
 White. **134800075**  
 3" Core, 40 yd roll  
 $\frac{1}{2}$ " Wide, each ..... \$9.80  
 Green. **134633005**    White. **134630005**  
 $\frac{3}{4}$ " Wide, each ..... \$13.00  
 Green. **134633075**    Yellow. **134632075**  
 $\frac{3}{4}$ " Wide, White, Pk/4. **134830075** ..... \$51.95  
 $\frac{3}{4}$ " Wide, Mixed, Pk/6. **134630600** ..... \$79.20  
 blue, green, orange, red, white & yellow



**Label Tape Dispenser, multi-roll. 134730000**..... \$41.35

**Test Strips** ..... *Call regarding availability*



**Fuses**

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. <b>312-100</b> .....	\$2.00
0.2A/250V. <b>312-200</b> .....	\$2.00
1A/250V. <b>312-001</b> .....	\$2.00
2A/250V. <b>312-002</b> .....	\$2.00
3A/250V. <b>312-003</b> .....	\$2.00
4A/250V. <b>312-004</b> .....	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. <b>313-700</b> .....	\$2.00
1A/250V. <b>313-001</b> .....	\$2.00
2.5A/250V .....	\$2.00

5 x 20 mm

Fast-Acting, GMA


0.5A/250V. <b>235-500</b> .....	\$2.00
1A/250V. <b>235-001</b> .....	\$2.00
3A/250V. <b>235-003</b> .....	\$2.00
5A/125V. <b>235-005</b> .....	\$2.00

Time-Delay (Slow-Blow), Glass tube

0.5A/250V. <b>239-500</b> .....	\$2.00
1A/250V. <b>239-001</b> .....	\$2.00
2A/250V. <b>239-002</b> .....	\$2.00
3A/250V. <b>239-003</b> .....	\$2.00
5A/250V. <b>239-005</b> .....	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. <b>215-001</b> .....	\$2.00
-------------------------------	--------



Cone Point



Cup Point

**Set Screws, steel**

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. <b>77702</b> .....	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. <b>77704</b> .....	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. <b>77706</b> .....	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. <b>77708</b> .....	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. <b>77722</b> .....	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. <b>77724</b> .....	\$0.50

**Tools**




**Cutter, diagonal, 5". PLR-766.00** ..... \$10.95



Demagnetize or magnetize parts and tools.

**Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00** ..... \$80.95

**Gauges**



Vernier Caliper measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

**Vernier Caliper, digital, 3-button. 600VCD001** ..... \$28.95



**Compasses, replaceable tips, Germany**

3" Long. <b>DIV-320.00</b> .....	\$11.95
6" Long. <b>DIV-326.00</b> .....	\$16.95
Tips (replacement), Pair. <b>DIV-320.05</b> .....	\$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

**Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A** ..... \$39.00



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

**Tubing Sizer, metric diameter gauge, PP. 133900000** ..... \$35.10