

Thermometers, Digital, Traceable®



Refrigerator



Refrigerator/
Freezer



| Model | Range | Resolution | Accuracy | Probe Length | Unit Size | Replacement Battery | Price |
|-------|--------------------------------|------------|----------|--------------|----------------------|---------------------|---------|
| 4352 | -58 to 302°F/ -50 to 150°C | 0.1° | ±0.2°C | 8" | 11 x 13/16" dia. | *357 | \$59.10 |
| 4354 | -58 to 572°F/ -50 to 300°C | 0.1° | ±0.5°C | 11 3/8" | 14 3/8 x 13/16" dia. | *357 | \$60.00 |
| 4371 | -58 to 572°F/ -50 to 300°C | 0.1° | ±1.0°C | 8" | 10 1/4 x 1 3/4" dia. | *LR1130 | \$28.65 |
| 4378 | -58 to 572°F/ -50 to 300°C | 0.1° | ±0.4°C | 8" | 10 1/4 x 1 3/4" dia. | * LR1130 | \$64.70 |
| 4158 | 14.2 to 131°F/ -9.9 to 55°C | 0.1° | ±0.5°C | N/A | 2 1/2 x 1 x 3/8" | * LR1130 | \$30.20 |
| 4127 | -58 to 158°F/ -50 to 70°C | 1° | ±1.0°C | 2 1/2" | 4 1/4 x 2 3/4 x 3/4" | *AA | \$53.90 |
| 4374 | -58 to 302°F/ -50 to 150°C | 0.1° | ±0.3°C | 5 1/4" | 6 1/2 x 2 1/8" dia. | *357 | \$53.25 |

*Battery (replacement): 357 / \$4.00; LR1130/ \$2.25. Warranty: 1 year



Thermometers, Glass, student grade

| Model | Range | Divisions | Length | Features | Price |
|---------|-------------------------------|-----------|--------|---|--------|
| 97-3104 | -20 to 230°F/ -30 to 110°C | 2°F/1°C | 5 3/4" | Red Alcohol/Aluminum V-Back. | \$2.20 |
| 97-3110 | -10 to 102°C | 1°C | 6" | Red Alcohol/White Background. Partial immersion. | \$3.25 |
| 97-3108 | 0-230°F/ -20 to 110°C | 2°F/1°C | 12" | Red Alcohol/White Background. Partial immersion. | \$3.25 |



Thermometer Rack, 25-place, polypropylene. 310THER01 \$26.25