

**Microscope Stage Laminate, USA. 16251** ..... \$8.00  
6.25 x 7.75" sheet provided with a #10 scalpel.  
White Teflon® (PTFE) on self-adhesive vinyl is 0.010" thick.



**Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420** ..... \$0.50



Mark slides, dishes, flasks, and chambers with a 1.8 mm ring for later referral. Works on upright and inverted microscopes.

**Object Marker, RMS thread (20.32 mm Dia.)**  
Black. **14361** ..... \$172.00  
Refill Ink, black, 10 mL, Shachihata. **22112** ..... \$7.00  
Adapters  
RMS/45 to CFI60, 15 mm long. **14111** ..... \$43.00  
RMS ID to 25 mm OD Thread, brass. **14115** ..... \$52.00

*Microscopes with Nikon CFI60 infinity-corrected ( $\infty$ ) objectives require a 14111 adapter.*  
*Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.*  
*Microscopes with 160 mm tube length objectives do not require an adapter.*



**Objectives, 160 mm tube length, DIN, Jenco, Meiji, National**  
Achromats, Motic/National  
4X, N.A. 0.10. **704-160** ..... \$32.00  
10X, N.A. 0.25. **710-160** ..... \$33.00  
40X, N.A. 0.65. **740-160** ..... \$42.00  
100X Oil, N.A. 1.25. **799-160** ..... \$69.00  
Super High Contrast Achromats, National  
4X, N.A. 0.10. **704-160ASC** ..... \$28.00  
10X, N.A. 0.25. **710-160ASC** ..... \$36.00  
40X, N.A. 0.65. **740-160ASC** ..... \$66.00  
100X Oil, N.A. 1.25. **799-160ASC** ..... \$99.00  
Semi-Plan Achromats, National (unless noted)  
2X, N.A. 0.05, Accu-Scope. **3193SP-00** ..... \$70.00  
4X, N.A. 0.10. **704-160SP** ..... \$32.00  
10X, N.A. 0.25. **710-160SP** ..... \$33.00  
40X, N.A. 0.65. **740-160SP** ..... \$42.00  
50X Oil, N.A. 0.95, Meiji. **MA694** ..... \$295.00  
60X, N.A. 0.85. **MA645** ..... \$315.00  
100X Oil, N.A. 1.25. **799-160SP** ..... \$69.00

Continued

Plan Achromats, National (unless noted)  
2.5X, N.A. 0.08, Meiji. **MA432** ..... \$285.00  
4X, N.A. 0.10. **704-160P** ..... \$33.00  
10X, N.A. 0.25. **710-160P** ..... \$55.00  
20X, N.A. 0.40. **720-160P** ..... \$125.00  
40X, N.A. 0.65. **740-160P** ..... \$104.00  
50X Oil, N.A. 0.95, Accu-Scope. **3190PL-00** ..... \$315.00  
60X, N.A. 0.85. **760-160P** ..... \$199.00  
100X Oil, N.A. 1.25. **799-160P** ..... \$168.00

Protect spare microscope objectives.



**Objective Case, DIN, polystyrene. 15601** ..... \$0.50




Extender tubes do not contain lenses.

**Objective Extender Tubes, RMS thread, brass, USA**  
10 mm Long. **14105** ..... \$44.20  
18 mm Long. **14106** ..... \$42.75

**Objective Shims, fit RMS thread (20.32 mm Dia.), USA**  
Parfocalize objectives with these stainless steel rings.  
25  $\mu$ m (0.001") Thick, Pk/10. **10032061** ..... \$16.15  
50  $\mu$ m (0.002") Thick, Pk/10. **10032062** ..... \$16.15  
75  $\mu$ m (0.003") Thick, Pk/10. **10032063** ..... \$16.15  
100  $\mu$ m (0.004") Thick, Pk/10. **10032064** ..... \$16.15  
125  $\mu$ m (0.005") Thick, Pk/10. **10032065** ..... \$16.15  
Assortment Pack, pk/10. **100320630** ..... \$19.80  
Includes two each of the following thicknesses: 25, 50, 75, 100, and 125  $\mu$ m.

Highly efficient, linear polarizing film has excellent color reproducibility. Transmission is 38% for a single polarizer, while extinction is better than 99.99% for crossed polars. Film is 0.76 mm thick.



**Polarizing Film, 50 x 50 mm Pair w/Plastic Case, USA. 15751** ... \$15.00