

# MICROSCOPES, CAMERAS & MAGNIFIERS

## Compound Microscopes

### Upright Microscopes



BK-211

#### Jenco

**BK-211** Upright Microscope, Binocular, Brightfield ..... \$1,215.00

- Widefield 10X/22 eyepieces
- Binocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- DIN, plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
- Ball-bearing mechanical stage is 19 cm(W) x 14 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- Centerable 1.25 N.A. Abbe condenser with rack and pinion focusing, and iris diaphragm
- 20 watt halogen Koehler-type illuminator
- 30.5 cm(D) x 20.3 cm(W) x 38.1 cm(H)
- 5-year mechanical & optical/2-year electrical warranty

**BK-311** Upright Microscope, Trinocular, Brightfield ..... \$1,357.00

- Trinocular version of BK-211



#### Meiji

**MT4200H** Upright Microscope, Binocular, Brightfield ..... \$2,439.00

- Widefield 10X/20 eyepieces
- Binocular, 30° inclined, Siedentopf design viewing head has single diopter; 53-75 mm interpupillary adjustment
- Infinity-corrected, plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
- Ball-bearing, ceramic coated mechanical stage is 17 cm(W) x 14 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- Centerable 1.25 N.A. Abbe condenser has rack and pinion focusing, and iris diaphragm
- 30 watt halogen Koehler-type illuminator
- 39 cm(D) x 21 cm(W) x 41 cm(H)
- Lifetime mechanical & optical/90-day electrical warranty

MT5300H



#### Meiji

**MT5200H** Upright Microscope, Binocular, Brightfield ..... \$2,695.00

- Widefield 10X/22 eyepieces
- Binocular, 30° inclined, Siedentopf design viewing head has single diopter; 53-75 mm interpupillary adjustment
- Infinity-corrected, plan semiapochromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
- Ball-bearing, ceramic coated mechanical stage is 19 cm(W) x 14 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- Centerable 1.25 N.A. Abbe condenser has rack and pinion focusing, and iris diaphragm
- 30 watt halogen Koehler-type illuminator
- 39 cm(D) x 21 cm(W) x 41 cm(H)
- Lifetime mechanical & optical/90-day electrical warranty

**MT5300H** Upright Microscope, Trinocular, Brightfield ..... \$2,977.00

- Trinocular version of MT5200H
- 38 mm ISO photoport includes 1X C-mount coupler
- 39 cm(D) x 21 cm(W) x 46.5 cm(H)



168-ASC



169-ASC

#### National

**168-ASC** Upright Microscope, Binocular, Brightfield ..... \$775.00

- Widefield 10X/18 eyepieces
- Binocular, 30° inclined, Jentsch design viewing head with dual diopters; 54-76 mm interpupillary adjustment
- DIN, super high contrast achromat 4X, 10X, 40XR, and 100XR objectives on 4-place, inward facing nosepiece
- Rackless x-axis mechanical stage is 14 cm(W) x 13.5 cm(D); up-stop adjustment is provided
- Coaxial coarse and fine focus controls; coarse has tension adjustment
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm, and swing-out filter holder
- 3 watt variable LED illuminator, corded
- 24 cm(D) x 18.5 cm(W) x 40 cm(H)
- Lifetime mechanical & optical/3-year electrical warranty

**169-ASC** Upright Microscope, Trinocular, Brightfield ..... \$989.00

- Trinocular version of 168-ASC
- 0.5X C-mount adapter included

*Continued*