

Spatulas, Micro, stainless steel

- Square End/Spoon End
 - 19 cm (7½"). **320-411** \$4.15
 - 23 cm (9"). **320-414** \$4.25
- Round End/Tapered End, 18 cm (7") long. **320-370** \$4.65
 - Round end: 8 x 51 mm; Tapered end: 6 x 51 mm
- V-Shaped End/Round End, 16 cm (6¼") long. **HS15906** \$5.45



Aluminum



Polycarbonate

Scoops

- Aluminum, ribbed handle, round
 - 350 mL (12 oz). **105ASC12** \$10.45
 - 700 mL (24 oz). **105ASC24** \$13.45
 - 1100 mL (38 oz). **105ASC38** \$17.95
 - 1700 mL (58 oz). **105ASC58** \$23.95
- Polycarbonate (PC), transparent, rectangular
 - 175 mL (6 oz). **1052882** \$9.95
 - 945 mL (32 oz). **1059F75** \$17.95
 - 1890 mL (64 oz). **1059F76** \$19.95



Stainless Steel

- V-Shaped (Reagent Digger), 16.5 cm (6½"). **600LSCSS6** \$1.75
- Handle, wooden, 13 cm long, each. **600LSHND1** \$2.05

Teasing Needles, Probes, Picks & Pin Vises



97-1302


Teasing Needles

- Angular, wooden handle, 14.5 cm (5¾"), Pk/12. **97-1302** \$4.30
- Straight, plastic handle, 14.5 cm (5¾"), Pk/12. **97-1313** \$4.05



- Stainless Steel, straight, 14 cm. **M28** \$5.00

Blunt tip is angled and tapered. Hexagonal handle is 6 mm wide.



Probes, Light Mall, 15.5 cm (6")

- Nickel Plated Steel. **15-800** \$2.90
(5800)
- Stainless Steel. **15-800SS** \$3.80



Probes, Bowman Lacrimal, double-ended


- Economy, sterling silver, 12.7 cm (5") long
 - Size 0000-000 (~0.4 – 0.5 mm). **24-706** \$11.65
 - Size 00-0 (~0.6 – 0.7 mm). **24-708** \$11.65
 - Size 1-2 (~0.9 – 1.0 mm). **24-710** \$11.65
 - Size 3-4 (~1.1 – 1.2 mm). **24-712** \$11.65



320-216A

Probe/Seekers, stainless steel

- Double-Ended, 15.5 cm. **320-216A** \$4.50
 - One end tapered to a point; one end curved with a rounded tip.
- Single End, 12.5 cm (5"). **320-217A** \$4.40
 - Curved, fine point tip with hexagonal handle.



Probe with Bulb Tips, 25 cm (10"), stainless, German. 541-807 \$17.50



Probe with Eye & Bulb Tip, flexible, 15.5 cm (6"). 23-496 \$2.25
(5850)

Picks (Explorers) *Call regarding availability*



PIN-220.00



PIN-212.00

Hold needles, probes, wires, drill bits, etc.

Pin Vise, Swivel-Head, 0-3.3 mm opening. PIN-220.00 \$7.95
Two double-ended collets are provided; one in screw chuck, one inside handle.

Pin Vise, Mini, 0-0.8 mm opening. PIN-212.00 \$4.95
blackened tool steel