



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA 5-1450 \$419.00
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$21.00



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.

Sterilizer, Infrared Heat, BactiZapper™. B1000 \$330.00





- Blue film/white paper
- EO gas/steam
- Latex-free
- Self-sealing

Sterilization Pouches


3½ x 9" (8.9 x 22.9 cm), Bx/200. **440166423** \$21.00
 5¼ x 10" (13.3 x 25.4 cm), Bx/200. **440166424** \$31.40
 7½ x 13" (19.1 x 33 cm), Bx/200. **440166425** \$55.50

6100A



Sterilizing Trays, plastic, autoclavable, PST, USA

B/L/M (base, lid & mat) included.
 6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A** \$23.00
 (500252)
 ★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A** \$33.00
 10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A** \$67.00
 (500254)




B//L/2M (base, insert tray, lid & 2 mats) included.
 15 x 10 x 1½" (38.1 x 25.4 x 3.8 cm). **1520A** \$405.00




Instrument Trays

Polypropylene, Deep w/Flat Cover, autoclavable
 8½ x 6 x 2¾" (22 x 15 x 7 cm). **60081736** \$10.10



Stainless Steel

Shallow (Mayo), **USA**
 13¾ x 9¾ x ¾" (34.8 x 24.9 x 1.9 cm). **13F** \$43.50
 15½ x 10½ x ¾" (38.4 x 27 x 1.9 cm). **15F** \$51.65
 Deep w/Dome Cover
 8¾ x 3¼ x 1½" (22.2 x 8.4 x 3.8 cm). **35-309** \$31.25
 Deep - Order trays and covers separately, **USA**
 8⅞ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **952** \$28.40
 Cover, Fits 952 Tray, dome w/handle. **952D-2** \$44.90
 17⅞ x 4½ x 2⅞" (43.6 x 11.2 x 5.4 cm). **170-0** \$50.75
 Cover, Fits 170-0 Tray, dome w/handle. **170-2** \$40.05



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x ½"

Sharpening Stone, fine grit (hard). ABR-800.10 \$23.95

Includes a plastic cleaning basket. Requires only 42 kHz of power. Dimensions: 6 x 3½ x 2"(D).



Ultrasonic Cleaner, 600 mL, mechanical timer. CLN-520.20 \$45.95
 Accessories

Cleaning Basket, fine mesh, 5 x 4 x 1¼". **CLN-650.10** \$14.95
 Cleaning Solution, 1 Quart. **CLN-850.10** \$9.95
 Biodegradable concentrate mixes 1:30.