

SLIDE & CASSETTE PROCESSING

Microscope Slides


Slides, general, 75 x 25 x 1 mm

	bx/72	cs/20 bx
Plain, ground edges		
Premier®, safety corners. 9101 	\$3.90	76.00
Globe, 90° corners, "Globe". 1301	—	101.00
Frosted, beveled edges		
Premier®, sandblasted, safety corners. 8205	\$6.05	115.00
Globe, etched, clipped corners, "Globe". 1334 ...	—	124.00
White End, clipped corners, "Premiere"		
Premier®. 9108W-E	—	115.00
Superfrost®, "SPECIMEN". 4951-001T	—	257.00
<i>(12-550-143, 48311-600)</i>		


Superfrost® Plus



Concavity, 2-Well





Slides, specialty


Charged (silanated), 75 x 25 x 1 mm		
Recommended for adhering tissue sections and cell preparations.		
Diamond™, white end, Globe		
Cs/20 Bx. 1358W		\$275.00
Superfrost® Plus, white end, Erie, USA		
Cs/20 Bx. 4951PLUS-001		\$597.00
Preferred slide for Ventana instrumentation.		
<i>(12-550-15, 48311-703)</i>		
Concavity		
75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy		
1-Well, Bx/12. 600CS3X11		\$6.45
2-Well, Bx/12. 600CS3X12		\$8.80
76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie		
1-Well, Bx/12. 1526		\$64.95
2-Well, Bx/12. 1527		\$72.55
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA		
QMS3X1X1 , each.....		\$20.00
Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA		
FSMS3X1X1 , each.....		\$29.00


Coverslips (Cover Glasses)



Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)		
Circles, Germany*		
25 mm, ~222 pcs/oz, Cs/10 oz. 0-25		\$214.00

#1 Thickness (0.13-0.17 mm)		
Squares & Rectangles, Erie, USA		
22 x 22 mm, ~162 pcs/oz. 1-2222	\$7.45	62.70
22 x 30 mm, ~119 pcs/oz. 1-2230	—	57.30
22 x 50 mm, ~71 pcs/oz. 1-2250	—	57.30
24 x 50 mm, ~65 pcs/oz. 1-2450	—	57.30
24 x 60 mm, ~54 pcs/oz. 1-2460	—	57.30

#1 Thickness (0.13-0.17 mm)		
Circles, Germany*		
12 mm, ~663 pcs/oz. 1-12	\$268.00	—
18 mm, ~296 pcs/oz. 1-18	\$150.00	—
25 mm, 100 pcs/pk, Cs/10 Pk. 1-25		\$212.00

#1.5 Thickness (0.16-0.19 mm). <i>Recommended for most users</i>		
Squares & Rectangles, Erie, USA		
18 x 18 mm, ~209 pcs/oz. 1.5-1818	—	61.40
22 x 22 mm, ~140 pcs/oz. 1.5-2222	\$7.45	62.70
22 x 40 mm, ~77 pcs/oz. 1.5-2240	—	57.30
22 x 50 mm, ~62 pcs/oz. 1.5-2250	—	57.30
24 x 30 mm, ~94 pcs/oz. 1.5-2430	—	57.30
24 x 40 mm, ~71 pcs/oz. 1.5-2440	—	57.30
24 x 50 mm, ~56 pcs/oz. 1.5-2450	\$6.85	57.30
24 x 60 mm, ~47 pcs/oz. 1.5-2460	—	57.30

#1.5 Thickness (0.16-0.19 mm). <i>Recommended for most users</i>		
Circles, Germany*		
12 mm, ~552 pcs/oz. 1.5-12	\$223.00	—
18 mm, ~245 pcs/oz. 1.5-18	\$250.00	—
22 mm, ~164 pcs/oz. 1.5-22	—	\$223.00
25 mm, ~127 pcs/oz. 1.5-25	—	\$252.00
Circles, Erie, USA		
40 mm, ~54 pcs/oz. 1.5-40	\$29.00	244.00
40 mm required for Biotech's FCS2 chamber.		

*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Counting Chamber - see Counting Chamber Coverslips in MICROSCOPE Accessories & Parts section

Fused Quartz, **USA**

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. QCS1X1	\$19.00
1" Circle, each. QCS1X0	\$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. R525025	\$32.40
25 mm Circle, each. R525000	\$23.40

Fused Silica, **USA**

0.15-0.25 mm Thick, Corning 7980

1" Square, each. FSCS1X1	\$21.50
1" Circle, each. FSCS1X0	\$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. R425025	\$32.00
25 mm Circle, each. R425000	\$28.80