

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker. Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot®, permanent, extra fine point	ea	bx/12
Black. 44102	\$1.70	17.45
Blue. 44103	\$1.70	17.45
Green. 44105	\$1.70	17.45
Red. 44104	—	17.45



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Pens, Super ^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic		
Regular Tip (4 mm tapered nib). K-600	\$61.00	
(71310, 195505)		
Fine Tip (2.5 mm tapered nib). K-650	\$46.00	
(71312, 195506)		

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

China Marker (Wax Pencil), Bx/12. 02089	\$11.50
--	---------



Diamond PerfectMark™, permanent, black, Bx/10. 1000	\$23.20
--	---------



Sharpie Retractable, permanent	ea	bx/12
Fine Tip		
Black. 36701	—	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip		
Black. 1735790	\$2.85	29.65
Blue. 1735792	—	29.65

Diamond-Tipped



Carbide-Tipped

Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$22.75
Diamond Tip (replacement). RDT-1	\$12.90




Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers

Small, 26 x 18 cm (10 x 7"). XH-2002	\$257.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$345.00

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

10-Slide, polystyrene, USA	ea	pk/3
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00



20-Slide, Simport		
Clear Lid (PETG), black base (ABS), each. M920-1	\$162.00	
Black Lid (ABS), black base (ABS), each. M920-2	\$162.00	

Continued



Incubators, MyTemp™ Mini Digital, Benchmark

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Two models: Heat Only, or Heat/Cool
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
14.5"(D) x 13.2"(W) x 18.7"(H) exterior
- 2-year warranty

Heat Only, 1° above ambient to 60°C. **H2200-H** \$426.00
 Heat/Cool, 15° below ambient to 60°C. **H2200-HC** \$656.00
 Heat/Cool/CO₂, 10° below ambient to 55°C. **H2300-HC2** \$3,182.00

Markers

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

China Marker (Wax Pencil), Bx/12. **02089** \$11.50



Sharpie Retractable, permanent			
Fine Tip		<u>ea</u>	<u>bx/12</u>
Black. 36701		---	\$29.65
Blue. 36703	\$2.85		29.65
Red. 36702	\$2.85		29.65
Ultra-Fine Tip		<u>ea</u>	<u>bx/12</u>
Black. 1735790	\$2.85		29.65
Blue. 1735792	---		29.65

Mortar and Pestle Sets

Glass, **USA**

59 mL (2 oz) Capacity. **826-2** \$15.75
 118 mL (4 oz) Capacity. **826-3** \$16.50
 237 mL (8 oz) Capacity. **826-4** \$19.65



Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. **LPC-154** \$4.70
 80 mL Capacity, 80 mm top OD. **LPC-154-1** \$6.30
 100 mL Capacity, 105 mm top OD. **LPC-155** \$7.75
 320 mL Capacity, 135 mm top OD. **LPC-156** \$10.55
 500 mL Capacity, 170 mm top OD. **LPC-157** \$15.30

Petri Dishes



Petri Dishes

- PS (polystyrene)
- 90 x 20 mm, economy
Pk/10, nonsterile. **51651** \$5.00
 - 35 x 10 mm, vented lid, stacking ring, Simport
Cs/500, sterile sleeve/20. **D210-15** \$133.00
 - 60 X 15 mm, stackable, Kord-Valmark
Cs/500, sterile sleeve/20. **2901** \$94.00
 - 100 X 15 mm, Kord-Valmark
Slippable (rimless lid for automation)
Standard, Cs/500, sterile sleeve/25. **2917** \$106.00
ISO Mark, Cs/500, sterile sleeve/25. **2910P** \$106.00
 - ★ Stackable
Standard, Cs/500, sterile sleeve/25. **2900** \$109.00
 - 100 X 25 mm, stackable, Kord-Valmark
Cs/340, sterile sleeve/17. **2906** \$109.00
 - Glass (borosilicate), reinforced rims, economy
90 x 15 mm, each. **31015003** \$2.80
100 x 15 mm, each. **31015004** \$3.20



Store (54) 100 mm petri dishes in this clear, 6-column rack. 34(L) x 20(W) x 23 cm(H).

Petri Dish Rack, 100 mm, acrylic. 189850100 \$126.00