


Tissue (Embedding) & Biopsy Cassettes




Base Molds, plastic, disposable, Simport

- 15 x 15 x 5 mm, Cs/1000. **M475-2** \$129.00
- 24 x 24 x 5 mm, Cs/1000. **M475-3** \$129.00
- ★ 30 x 24 x 5 mm, Cs/1000. **M475-4** \$129.00
- 37 x 24 x 5 mm, Cs/1000. **M475-5** \$129.00
- 37 x 24 x 10 mm, Cs/500. **M475-10** \$67.30

Designed for Simport Macrosette™.

Solvent resistant,
0.2 x 0.2 mm mesh.



Nylon




Paper

Biopsy Bags

- Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011** \$235.00
- Paper (Tea Bags), 2 1/4 x 2 1/4", Pk/500. **56100** \$38.00

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.




M476-2

Biopsy Pads

- Rectangular
 - Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014** \$40.00
 - Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1** \$190.00
- Circular, Simport
 - Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2** \$231.00

Biopsy Cassettes




Economy Biopsy Cassettes

Biopsy Cassettes, Detached Lids (packed separately)

Economy, White, small square pores, Pk/500. **BC-0109W** \$42.90



Simport Histosette I



Simport Micromesh

Biopsy Cassettes, Attached Lids, acetyl polymer

Histosette™ I, 1 mm square pores, Pk/1500 1 pk pk/3
\$126.00 355.00
30° front labeling area.


- Blue. **M491-6** Pink. **M491-3**
- Green. **M491-4** White. **M491-2**

Uniset™, hinged lids, Pk/1500 1 pk pk/3
\$135.00 383.00
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

- Blue. **M506-6** Pink. **M506-3**
- Green. **M506-4** White. **M506-2**

Micromesh™, 0.38 mm pores, Pk/1000 1 pk pk/3
\$211.00 589.00
45° front labeling area.

- Blue. **M507-6** Pink. **M507-3** Yellow. **M507-5**
- Green. **M507-4** White. **M507-2**



Compatible with Primera Signature™ and Tissue-Tek SmartWrite™ cassette printers.

Biopsy Cassettes, Taped

Micromesh™, 0.38 mm pores, Pk/2000 \$353.00
45° front labeling area.

- Blue. **M407-6T** Pink. **M407-3T** Yellow. **M407-5T**
- Green. **M407-4T** White. **M407-2T**

Continued

