



Bottle Top Dispensers

- 0.5-5 mL. **600BTDR2**\$251.00
- 5-60 mL. **600BTDR5**\$298.00

Boxes



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white
 2 1/8 x 1 1/8 x 5/8", Pk/72. **20047**\$31.10

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**\$121.00
 Each cardboard box can store approximately 540 paraffin blocks.



246530002

246530001

Boxes, Glass Disposal, cardboard, polyethylene bag included
 Bench Model, 8 x 8 x 10", Pk/6. **246530002**\$67.55
 Floor Model, 12 x 12 x 27", Pk/6. **246530001**\$93.80

Hinged
 Polystyrene




Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80


BRS-969.10



LB-40



LB-38




LB-30



LB-14



Brushes

- Bench, horse hair, 10" long. **BRS-969.10**\$9.20 ↓
 Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".
- Brass Wire, 9 3/4" long. **LB-40**\$8.70
 Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".
- Flat Bristles, paint style
 Camel hair
 1/2" Wide, 6 1/4" total length. **LB-36**\$3.90
 1" Wide, 7" total length. **LB-37**\$6.15
 2" Wide, 7 3/4" total length. **LB-38**\$5.30
-  Horsehair, heavy-duty bristles, metal handle. **33761880**\$0.85
- Round Bristles, camel hair, 8 3/4" long. **LB-30**\$10.25
- Swivel-Head, black nylon bristles. **LB-14**\$4.20
 Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



- Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany
 Brass. **BRS-295.00**\$12.20
 Nylon. **BRS-293.00**\$9.70
 Steel. **BRS-292.00**\$12.80

Burners



- Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313**\$6.80
 Wicks (replacement), Pk/12. **97-5310**\$5.00



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



- Burners, Bunsen**, natural gas, Humboldt Mfg, **USA**
 Standard, flame stabilizer. **H-5870**\$22.45
 High Temperature, Meeker-type. **H-5550**\$78.90

Continued