

## Student Microscopes



### National

**185 Student Microscope, Monocular** ..... \$89.00

- Upright and unreversed image
- Widefield 10X eyepiece
- Monocular, 45° inclined viewing head rotates 360°
- 2X objective provides 20X magnification
- Stage has two clips and reversible black/white plate
- Rack and pinion focusing
- 26 cm (10¼") tall
- 5-year warranty
- Grades K-12



### C&A Scientific

**MFL-05 Student Microscope, Monocular** ..... \$130.00

- Widefield 10X eyepiece with pointer
- Monocular, 45° inclined viewing head rotates 360°
- Achromat 4X, 10X and 40XR objectives
- Stage provided with safety stop to protect slides & objectives
- Coarse focus knob
- In-stage condenser with 6-place disc diaphragm
- Adjustable intensity LED illuminator, equivalent to a 20 watt tungsten bulb, operates on rechargeable battery or power cord
- 27.8 cm (11") tall
- 1-year warranty
- Five prepared glass slides
- Grades K-8



### National

**131-CLED Student Microscope, Monocular, DIN**..... \$237.00

- Widefield 10X eyepiece w/ pointer accepts optional reticle
- Monocular, 45° inclined viewing head rotates 360°
- DIN achromat 4X, 10X, and 40XR objectives
- Stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- LED illuminator, corded (equivalent to a 20 watt tungsten bulb)
- 36 cm (14") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K-12



### National

**132-CLED-MS Student Microscope, Dual View, DIN**..... \$378.00

- Widefield 10X eyepieces w/ pointers accept optional reticle
- Vertical & 30° inclined tubes mounted on 360° rotating head
- DIN achromat 4X, 10X, and 40XR objectives
- Graduated mechanical stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- LED illuminator, corded (equivalent to a 20 watt tungsten bulb)
- 40 cm (15¼") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K-12