



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

- Bag Sealer**, electric, 12" (30.5 cm) Wide. **HS28272**..... \$183.00
- Bag Sealer Replacement Part Kit. **HS28280**..... \$31.00
- Poly Tubing
  - 4 mil
    - 4" x 1500' Roll. **1143** ..... \$110.00
    - 6" x 1500' Roll. **1144** ..... \$162.00
  - 6 mil
    - 4" x 1000' Roll. **1658** ..... \$132.00
    - 6" x 1000' Roll. **1659** ..... \$162.00
- Poly Tubing Dispenser. **H-28**..... \$34.00

### Balances

#### Balance, Portable



|                        |                                 |
|------------------------|---------------------------------|
| <b>Model</b>           | <b>DWP-1001</b>                 |
| <b>Capacity (g)</b>    | 1000                            |
| <b>Readability (g)</b> | 0.1                             |
| <b>Modes</b>           | g, oz, ct, gn                   |
| <b>Pan</b>             | 14.6 cm (5.75") stainless steel |
| <b>Power</b>           | 4 AA/AC adapter                 |
| <b>Price</b>           | \$72.00                         |

Warranty: 1 year

#### Balances, Topping

Benchmark  
W3200-1200



| Model                  | W3200-500         | W3200-1200   |
|------------------------|-------------------|--------------|
| <b>Capacity (g)</b>    | 500               | 1200         |
| <b>Readability (g)</b> | 0.001             | 0.01         |
| <b>Modes</b>           | g, oz, ct         | g, oz, ct    |
| <b>Pan</b>             | 11.5 cm (4.5")    | 16 cm (6.3") |
| <b>Draft Shield</b>    | Yes, 3-door glass | No           |
| <b>Price</b>           | \$1,040.00        | \$842.00     |

Warranty: 2 years

### Balances, Analytical

Benchmark



| Model                  | ★ W3100-120                         | W3100A-120                        |
|------------------------|-------------------------------------|-----------------------------------|
| <b>Calibration</b>     | External<br>(100 g weight provided) | Internal<br>(Quick-Cal™ software) |
| <b>Capacity (g)</b>    | 120                                 | 120                               |
| <b>Readability (g)</b> | 0.0001                              | 0.0001                            |
| <b>Modes</b>           | g, oz, lb, ct                       | g, oz, lb, ct                     |
| <b>Pan</b>             | 9 cm (3.5")                         | 9 cm (3.5")                       |
| <b>Price</b>           | \$1,548.00                          | \$1,645.00                        |

Warranty: 2 years

### Baths



**Dry Baths**, digital, Benchmark (2-year warranty)

- 1-Block Capacity. **BSH1001**..... \$372.00
- 2-Block Capacity. **BSH1002**..... \$446.00
- 4-Block Capacity. **BSH1004**..... \$753.00

*Blocks sold separately – please specify configuration*

Blocks, aluminum

- Holds (48) 0.2 mL Tubes. **BSW02**..... \$93.00
- Holds (24) 0.5 mL Tubes. **BSW05**..... \$93.00
- Holds (24) 1.5 mL Tubes. **BSW1500**..... \$93.00
- Holds (24) 1.5/2.0 mL Tubes. **BSW1520**..... \$93.00
- Holds (20) 12/13 mL Tubes. **BSW13**..... \$93.00
- Holds (12) 15 mL Tubes. **BSW15**..... \$93.00
- Holds (5) 50 mL Tubes. **BSW50**..... \$93.00



Racks are only included with model B2000-4.

**Water Baths**, digital, MyBath™, Benchmark (2-year warranty)

- 4 L, 3 tube racks. **B2000-4** ..... \$769.00
- 8 L. **B2000-8**..... \$832.00
- 12 L. **B2000-12**..... \$995.00