



Non-Vacuum

Glazed porcelain sample plate has 3.5-4.5 cm center hole.

Vacuum



Desiccators

- Non-Vacuum, glass, economy
24 cm (10") ID. **125720002**.....\$65.75
- Vacuum, polycarbonate top/polypropylene bottom
138 mm Plate Diameter. **60055204**.....\$53.80
- 190 mm Plate Diameter. **60055205**.....\$73.25
- 238 mm Plate Diameter. **60055206**.....\$92.60



- Dippers**, 500 mL polypropylene cup with HDPE handle
- 3' Handle. **367840016**\$53.75
 - 6' Handle. **367850016**\$70.00
 - Replacement Cup, screw-on. **367860016**\$12.80



Autoclavable polypropylene bottle with square bottom holds 170 mL (6 oz). Simport

- Drosophila Stock Bottle**, Cs/500. **B352**.....\$144.00



188050000



HS23243A



188190023

Drying Racks

- Non-Electric
Bel-Art. **188050000**.....\$105.00
Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 17 x 43 x 48 cm(H).
- Heathrow. **HS23243A**.....\$115.00
- Electric, double-sided, 3 tier, Bel-Art. **188190023**.....\$1,045.00
1-year warranty
- Dust Filter & Holder (replacement). **999321563**.....\$12.00
- Heating Element (replacement). **999321564**.....\$193.00

Eyewear

Absorbs 99.9% UV



- Glasses**, Safety, integral side shield. **43574002**.....\$5.60



Flexible frame and adjustable elastic strap provides a tight seal.

- Goggles**, Safety, vented (chemical splash). **97-2001**.....\$3.70

Film & Foil



- Parafilm**, thermoplastic sealing film, **USA**
- 2" x 250' (5 cm x 76 m) Roll. **HS234526A**.....\$27.00
 - 4" x 125' (10 cm x 38 m) Roll. **HS234526B**.....\$27.50
 - 4" x 250' (10 cm x 76 m) Roll. **HS234526C**.....\$45.95

Holds one or two 2" rolls, or one 4" roll.



- Parafilm Dispenser**, ABS plastic, blue. **HS234525B**.....\$39.35



- Foil**, **Super-Polyfoil™**, polyester-backed aluminum foil
- 4 x 4" (10.2 x 10.2 cm) Sheets, Pk/100. **132000044**.....\$35.70
 - 5 1/8 x 5 1/8" (13 x 13 cm) Sheets, Pk/100. **132000066**.....\$46.55