

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610 (Whatman 5; S&S 602)	
9 cm, Pk/100. 6100-0900	\$9.70
12.5 cm, Pk/100. 6100-1250	\$12.25
15 cm, Pk/100. 6100-1500	\$14.75
20.5 cm, Pk/100. 6100-2050	\$23.30
24 cm, Pk/100. 6100-2400	\$33.05
Medium flow (30 mL/min), 2 µm retention, Grade 642 (Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)	
9 cm, Pk/100. 6420-0900	\$7.55
15 cm, Pk/100. 6420-1500	\$13.70
Medium flow (55 mL/min), 2.5 µm retention, Grade 601 (Whatman 1, 1CHR; S&S 593, 597, 2043A)	
4.25 cm, Pk/100. 6010-0425	\$5.80
9 cm, Pk/100. 6010-0900	\$7.80
12.5 cm, Pk/100. 6010-1250	\$11.30
15 cm, Pk/100. 6010-1500	\$13.75
18.5 cm, Pk/100. 6010-1850	\$19.95
20.5 cm, Pk/100. 6010-2050	\$21.50
33 cm, Pk/100. 6010-3300	\$58.05
Medium flow (60 mL/min), 6 µm retention, Grade 613 (Whatman Student; S&S 595, 0903)	
9 cm, Pk/100. 6130-0900	\$6.65
15 cm, Pk/100. 6130-1500	\$11.65
20.5 cm, Pk/100. 6130-2050	\$17.50
33 cm, Pk/100. 6130-3300	\$38.90
Medium flow (200 mL/min), 10 µm retention, Grade 631 (Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)	
20.5 cm, Pk/100. 6310-2050	\$21.70
33 cm, Pk/100. 6310-3300	\$58.60
Fast flow (235 mL/min), 25 µm retention, creped, Grade 615	
9 cm, Pk/100. 6150-0900	\$6.85
12.5 cm, Pk/100. 6150-1250	\$10.00
15 cm, Pk/100. 6150-1500	\$12.00
20.5 cm, Pk/100. 6150-2050	\$18.00
33 cm, Pk/100. 6150-3300	\$40.00
Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617	
9 cm, Pk/50. 6170-0900	\$6.25
12.5 cm, Pk/50. 6170-1250	\$8.40
15 cm, Pk/50. 6170-1500	\$9.80
20.5 cm, Pk/50. 6170-2050	\$14.90
33 cm, Pk/50. 6170-3300	\$31.90

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542 (Whatman 2V)	
12.5 cm, Pk/100. 5420-1250	\$65.15
18.5 cm, Pk/100. 5420-1850	\$73.75



Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, Grade 909 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)	
32 cm, Pk/100. 9090-3200	\$63.25
Fast flow (245 mL/min), 27 µm retention, creped, Grade 962 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)	
11 cm, Pk/100. 9620-1100	\$9.40

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, Grade 94 (Whatman 42; S&S 589/3)	
9 cm, Pk/100. 0940-0900	\$27.15
15 cm, Pk/100. 0940-1500	\$46.70
Medium flow (30 mL/min), 2 µm retention, Grade 74 (Whatman 40; S&S 589/2, 589/4)	
9 cm, Pk/100. 0740-0900	\$27.15
15 cm, Pk/100. 0740-1500	\$46.70
Fast flow (220 mL/min), 15 µm retention, Grade 55 (Whatman 54, 541; S&S 1505)	
9 cm, Pk/100. 0550-0900	\$30.65
15 cm, Pk/100. 0550-1500	\$49.25

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, Grade 601 (Whatman 1, 1CHR; S&S 593, 597, 2043A)	
46 x 57 cm, Pk/100. 6018-4657	\$110.00
0.34 mm thick, 140 mm/30 min absorption, Grade 238 (Whatman 3MM, S&S 598, 2316)	
7 x 10 cm, Pk/100. 2388-0710	\$26.05
20 x 20 cm, Pk/100. 2388-2020	\$57.25
46 x 57 cm, Pk/100. 2388-4657	\$209.00
0.42 mm thick, 150 mm/30 min absorption, Grade 237 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703)	
46 x 57 cm, Pk/100. 2378-4657	\$240.00
0.95 mm thick, 180 mm/30 min absorption, Grade 243 (Whatman 17CHR; S&S 2668)	
11 x 14 cm, Pk/100. 2438-1114	\$24.75
20 x 20 cm, Pk/100. 2438-2020	\$45.95
46 x 57 cm, Pk/100. 2438-4657	\$352.00
2.48 mm thick, 240 mm/30 min absorption, Grade 320 (S&S 2727)	
11 x 14 cm, Pk/50. 3208-1114	\$39.15
16 x 18 cm, Pk/50. 3208-1618	\$54.55
20 x 20 cm, Pk/50. 3208-2020	\$66.55



Smooth tips &
polypropylene grips

Forceps, Filter, 11.5 cm. 300-025	\$4.40
--	--------