


C-Shaped



O-Shaped

**Flask Weights**, vinyl-coated lead

C-Shaped lead rings, blue, Heathrow

Fits 250-1000 mL Flasks, 51 mm ID, 472 g. **HS8881B** ..... \$28.90Fits 500-2000 mL Flasks, 54 mm ID, 599 g. **HS8881C** ..... \$31.85

O-Shaped lead ring, blue, Heathrow

Fits 1000-4000 mL Flasks, 70 mm ID, 970 g. **HS8882D** ... \$32.00Lead Wire, 300 mm long, 114 g. **183240000** ..... \$32.65**Freezer Racks**

175HS-9-3



175HD-442



175V-13-2L

**Freezer Racks**, stainless steel, **USA****Horizontal** (fits upright freezers)Horizontal, Holds 2" Boxes 1-9    ≥109-Place (3L x 3H). **175HS-9-2** ..... \$54.00    48.00

Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)

★ 12-Place (3L x 4H). **175HS-12-2** ..... \$72.00    64.00

Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)

15-Place (3L x 5H). **175HS-15-2** ..... \$90.00    80.00

Dimensions: 16.31(L) x 5.5(W) x 11"(H)

16-Place (4L x 4H). **175HL-16-2** ..... \$96.00    85.00

Dimensions: 22(L) x 5.5(W) x 8.75"(H)

20-Place (4L x 5H). **175HL-20-2** ..... \$122.00    108.00

Dimensions: 22(L) x 5.5(W) x 11"(H)

Horizontal, Holds 3" Boxes, fits upright freezers

9-Place (3L x 3H). **175HS-9-3** ..... \$56.00    50.00

Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)

12-Place (4L x 3H). **175HL-12-3** ..... \$74.00    66.00

Dimensions: 22(L) x 5.5(W) x 9.5"(H)

**Horizontal Drawer** (fits upright freezers)

Horizontal Drawer, stainless steel drawers

Holds 2" Boxes 1-9    ≥109-Place (3L x 3H). **175HD-332** ..... \$79.00    70.00

Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)

12-Place (3L x 4H). **175HD-342** ..... \$103.00    92.00

Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)

15-Place (3L x 5H). **175HD-352** ..... \$132.00    117.00

Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)

16-Place (4L x 4H). **175HD-442** ..... \$141.00    125.00

Dimensions: 22(L) x 5.5(W) x 9.5"(H)

20-Place (5L x 4H). **175HD-542** ..... \$175.00    156.00

Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)

Continued

## Horizontal Drawer (cont)

Holds 5" Boxes 1-9    ≥106-Place (3L x 2H). **175RD-325** ..... \$69.00    61.00

Dimensions: 18.125(L) x 6(W) x 10.75"(H)

8-Place (4L x 2H). **175HD-425** ..... \$92.00    82.00

Dimensions: 24(L) x 6(W) x 10.75"(H)

Holds 15 mL tubes

160-Place (80L x 2H). **175HD-32ML-15** ..... \$93.00    82.00

Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)

220-Place (110L x 2H). **175HD-42ML-15**... \$120.00    107.00

Dimensions: 22(L) x 5.5(W) x 9.5"(H)

Holds 50 mL tubes

84-Place (42L x 2H). **175HD-42ML-50** ..... \$120.00    107.00

Dimensions: 22(L) x 5.5(W) x 9.5"(H)

**Horizontal Drawer-Tray** (fits upright freezers)

Drawer-Tray, aluminum drawer-trays

Holds 2" Boxes 1-9    ≥1012-Place (3L x 4H). **175HDT-342** ..... \$89.00    79.00

Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)

18-Place (3L x 6H). **175HDT-362** ..... \$133.00    118.00

Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)

24-Place (4L x 6H). **175HDT-462** ..... \$174.00    154.00


Dimensions: 22(L) x 5.5(W) x 12.75"(H)

**Vertical** (fits chest freezers)

Vertical, Holds 2" Boxes

13-Place. **175V-13-2L** ..... \$81.00    72.00

Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)

**Tube Storage Boxes** (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section**Funnels****Funnels**

Disposable, polypropylene, short stems, Simport

★ 55 mm Top ID, takes 11 cm filter, Cs/100. **F490-3** ..... \$31.75★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. **F490-4** ..... \$31.75

Glass (borosilicate), long stem (150 mm)

50 mm Top ID, each. **600GF614050** ..... \$6.3065 mm Top ID, each. **600GF614065** ..... \$7.6075 mm Top ID, each. **600GF614075** ..... \$7.80

Continued

Heavy-Duty, Molded-In Handle, polypropylene 250 mm Top ID, 2.5 L capacity, each. <b>60057232</b> .....	\$20.00
Micro, polypropylene, Bel-Art 24 mm Top ID, 2.6 mL capacity, Pk/12. <b>1146850024</b> .....	\$31.25
28 mm Top ID, 5.2 mL capacity, Pk/12. <b>146850028</b> .....	\$25.60
35 mm Top ID, 8.1 mL capacity, Pk/12. <b>146850035</b> .....	\$23.60
43 mm Top ID, 15.7 mL capacity, Pk/12. <b>146850043</b> .....	\$54.50
Powder, polypropylene Parallel Stem 80 mm Top ID, Pk/12. <b>146600080</b> .....	\$24.15
100 mm Top ID, Pk/6. <b>146600100</b> .....	\$23.50
145 mm Top ID, Pk/6. <b>60057224</b> .....	\$18.90
Tapered Stem 101 mm Top ID, each. <b>146810000</b> .....	\$29.65
152 mm Top ID, each. <b>146820000</b> .....	\$33.65
Standard Stem, Premium, polypropylene, Bel-Art Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock. 55 mm Top ID, takes 9 cm filter, Pk/12. <b>146760055</b> .....	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. <b>146760065</b> .....	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. <b>146760075</b> .....	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. <b>146760100</b> .....	\$32.60
Standard Stem, Economy, polypropylene 50 mm Top ID, each. <b>310FUNPP50</b> .....	\$0.80
62 mm Top ID, each. <b>310FUNPP62</b> .....	\$0.95
75 mm Top ID, each. <b>310FUNPP75</b> .....	\$1.15
100 mm Top ID, each. <b>310FUNPP100</b> .....	\$1.65

These 4 x 4" nonsterile  
gauze sponges are 12-ply.



**Gauze Pads**, Cs/2000, 10 pk/200. **44122000** .....

Sani-Cloth Plus



**Germicidal Wipes**, Tub/160. **440Q89072** .....

**Gloves & Hand Protection**

Cotton Gloves



Heat-Resistant Gloves



Nitrile Gloves



**Gloves**

Cotton

Protect optics and photographs from fingerprints.  
Can also be used as glove liners.

Small, Pk/12. **GLV-190.10** .....

Large, Pk/12. **GLV-190.20** .....

Heat-Resistant

Terry cloth gloves protect hands from high temperatures.

33 cm (13") Long, Pair. **132010000** .....

47 cm (18.5") Long, Pair. **132010001** .....

Latex, powder-free, 5 mil, fully textured **1 bx** **cs/10 bx**

X-Small, Bx/100. **73101** .....

Small, Bx/100. **73103** .....

Medium, Bx/100. **73105** .....

Large, Bx/100. **73107** .....

X-Large, Bx/100. **73109** .....

XX-Large, Bx/100. **73111** .....

Nitrile, powder-free, 5.9 mil, textured fingers, Halyard

Nitrile gloves are hypoallergenic and provide greater puncture  
and chemical resistance than latex or vinyl gloves.



**1 bx** **cs/10 bx**

X-Small, Bx/100. **73121** .....

Small, Bx/100. **73123** .....

Medium, Bx/100. **73125** .....

Large, Bx/100. **73127** .....

X-Large, Bx/90. **73129** .....

XX-Large, Bx/90. **73131** .....



**Glove Box Holder**, clear acrylic, Heathrow, **USA**

1-Box, 29 x 11 x 13 cm(H). **HS23471A** .....

2-Box, 29 x 11 x 26 cm(H). **HS23471B** .....



Safely grip hot or cold labware  
with this silicone rubber hand  
protector partially lined with  
mini suction cups. Withstands  
temperatures from 260°C (500°F)  
to -57°C (-70°F).

**Hand Protector**, red, Hot Hand@. **380000000** .....

**Hot Plate, Magnetic Stirrer & Hot Plate Stirrers**

1257225



H4000-H

**Hot Plates**, analog

6" Diameter Cast Iron Burner, 1000W. **1257225** .....

7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H** .....

Continued

Continued