




- Ceramic platform is 7.1 x 7.1" (18 x 18 cm)
- Stir speed is adjustable from 60-1500 rpm
- 2-year warranty

Magnetic Stirrer, analog, Benchmark. **H4000-S** \$219.00

Hot Plate Stirrer, 80-380°C, Benchmark. **H4000-HS** \$328.00

Homogenizers & Tissue Grinder



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizer, Handheld. D1000 \$937.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.
Benchmark BeadBlaster™ 24

Homogenizer, Benchtop, Microtube. D2400 \$6,796.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

0.1 mm, Pk/50. **D1031-01** \$101.00

0.5 mm, Pk/50. **D1031-05** \$101.00


Zirconium Beads, acid washed/heat treated/nuclease-free

1.5 mm, Pk/50. **D1032-15** \$141.00

3.0 mm, Pk/50. **D1032-30** \$141.00

Stainless Steel Beads, acid washed

2.8 mm, Pk/50. **D1033-28** \$141.00




Tissue Grinder, Dual® Glass Pestle

1 mL Capacity, **LG10640100** \$82.50

3 mL Capacity, **LG10640102** \$92.85

5 mL Capacity, **LG10640104** \$105.00



Hydrometers Access www.belart.com
and call us to place an order

Ice Buckets & Pans



Temperature Range: -196 to 100°C

Ice Buckets with lids, expanded PVC foam, Bel-Art, **USA**

2.5 L Capacity (Handy)

Black. **188482002** \$54.55 ↓

Blue. **188482001** \$54.55 ↓

Purple. **188482006** \$54.55 ↓

4 L Capacity (N'icer)

Black. **188484002** \$66.50 ↓

Blue. **188484001** \$66.50 ↓

Ice Pans w/o lids, expanded PVC foam, Bel-Art, **USA**

1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)

Black. **188481102** \$40.65 ↓

Blue. **188481101** \$40.65 ↓

Purple. **188481106** \$40.65 ↓

4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)

Black. **188484102** \$64.80 ↓

Blue. **188484101** \$64.80 ↓


Purple. **188484106** \$64.80 ↓

9 L Capacity (Maxi), 44 x 35 x 13 cm(H)

Black. **188489102** \$76.20 ↓

Blue. **188489101** \$76.20 ↓

Purple. **188489106** \$76.20 ↓



Tube Racks, Floating, PVC Foam, Bel-Art, USA

Circular. *Recommended for ice buckets*

44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$19.45 ↓

Black/**188480012** Green/**188480014** Red/**188480013**

30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$20.40 ↓

Black/**188480022** Green/**188480024** Red/**188480023**

Rectangular. *Recommended for ice pans*

36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm \$19.45 ↓

Black/**188480112** Green/**188480114** Red/**188480113**

24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm \$20.40 ↓

Black/**188480122** Green/**188480124** Red/**188480123**