

Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®
 40 mm (1½") Long. **371160002** \$25.70
 50 mm (2") Long. **371160003** \$38.40

Micro, Spinbar®

White

5 x 2 mm. **371190005** \$5.10
 7 x 2 mm. **371190007** \$5.10
 10 x 3 mm. **371190010** \$5.65
 12.7 x 3 mm. **371190127** \$5.65
 15 x 1.5 mm. **371190015** \$5.00

Red

7 x 2 mm. **371210013** \$5.65
 8 x 1.5 mm. **371210001** \$5.65

Octagonal, molded-on pivot ring, Spinbar®

White

12.7 x 8 mm (½ x 5/16"). **371101256** \$4.80
 15.9 x 8 mm (5/8 x 5/16"). **371100058** \$4.80
 25.4 x 9.5 mm (1 x 3/8"). **371100138** \$6.35
 38.1 x 9.5 mm (1½ x 3/8"). **371101128** \$7.50
 50.8 x 9.5 mm (2 x 3/8"). **371100238** \$10.00
 63.5 x 9.5 mm (2½ x 3/8"). **371102128** \$11.55
 76.2 x 12.7 mm (3 x ½"). **371100003** \$20.75

Red, 12.7 x 3.2 mm (½ x 1/8"), no pivot ring. **371090034** \$5.85

Stirring Bar Assortment, octagonal, Pk/6. **371600000** \$37.00

Includes two of each size: 25.4 x 8 mm (1 x 5/16"),
 38.1 x 8 mm (1½ x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

12 x 4.5 mm (7/16 x 3/16"). **371220012** \$5.05 ↓

Triangular, vanes w/cross-mounted magnets, Spinvane®


Fits 1 mL Vials. **371340000** \$9.95
 Fits 3/5 mL Vials. **371350000** \$9.95



Stirring Bar Rack, Spinbar® Garage. **377900000** \$48.95




Stirring Bar Retainer, Spinbar® Restrainer. **377861000** \$13.75




Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.15
 45.7 cm (18") Long. **377730000** \$17.55



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. **377870000** \$29.55



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.35
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$5.95
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.20
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$23.25