

Stoppers

Stoppers, rubber, black



Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.10	13.25
Size 3, ~29 pcs/lb. 50030	\$4.15	13.25
Size 4, ~23 pcs/lb. 50040	\$5.25	13.25
Size 5, ~20 pcs/lb. 50050	—	13.25
Size 6, ~15 pcs/lb. 50060	—	13.25
Size 6.5, ~14 pcs/lb. 500650	—	13.25
Size 7, ~11 pcs/lb. 50070	—	13.25
Size 8, ~9 pcs/lb. 50080	—	13.25
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.35	13.25
Size 2, ~43 pcs/lb. 50021	\$2.80	13.25
★ Size 3, ~30 pcs/lb. 50031	\$4.00	13.25
Size 4, ~24 pcs/lb. 50041	\$5.00	13.25
Size 5, ~21 pcs/lb. 50051	—	13.25
Size 6, ~18 pcs/lb. 50061	—	13.25
Size 9, ~8 pcs/lb. 50091	—	13.25
Size 12, ~3 pcs/lb. 50121	—	13.25
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$13.25	
Size 7, ~12 pcs/lb. 50072	\$13.25	
Size 8, ~9 pcs/lb. 50082	\$13.25	
Size 9, ~7 pcs/lb. 50092	\$13.25	
Size 10, ~5 pcs/lb. 50102	\$13.25	
Assortment, 1 lb. 50012	\$17.25	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$13.50

Non-fenestrated towels are sterile, latex-free and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515 \$25.50

Tape

Tape



Autoclave, steam


¾" x 60 Yd Roll. **44073ST048** \$8.50
 1" x 60 Yd Roll. **44073ST036** \$10.35↓

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll
 ½" Wide, each \$5.50
 Green. **134623005** Yellow. **134622005**
 ¾" Wide, Pk/4 \$23.75
 Green. **134820075** Yellow. **134810075**
 White. **134800075**
 3" Core, 40 yd roll
 ½" Wide, each \$10.05
 Green. **134633005** White. **134630005**
 ¾" Wide, each \$13.35
 Green. **134633075** Yellow. **134632075**
 ¾" Wide, White, Pk/4. **134830075** \$53.35
 ¾" Wide, Mixed, Pk/6. **134630600** \$81.20
 blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000 \$42.20