

Stoppers

Stoppers, rubber, black



Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.10	13.25
Size 3, ~29 pcs/lb. 50030	\$4.15	13.25
Size 4, ~23 pcs/lb. 50040	\$5.25	13.25
Size 5, ~20 pcs/lb. 50050	—	13.25
Size 6, ~15 pcs/lb. 50060	—	13.25
Size 6.5, ~14 pcs/lb. 500650	—	13.25
Size 7, ~11 pcs/lb. 50070	—	13.25
Size 8, ~9 pcs/lb. 50080	—	13.25
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.35	13.25
Size 2, ~43 pcs/lb. 50021	\$2.80	13.25
★ Size 3, ~30 pcs/lb. 50031	\$4.00	13.25
Size 4, ~24 pcs/lb. 50041	\$5.00	13.25
Size 5, ~21 pcs/lb. 50051	—	13.25
Size 6, ~18 pcs/lb. 50061	—	13.25
Size 9, ~8 pcs/lb. 50091	—	13.25
Size 12, ~3 pcs/lb. 50121	—	13.25
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$13.25	
Size 7, ~12 pcs/lb. 50072	\$13.25	
Size 8, ~9 pcs/lb. 50082	\$13.25	
Size 9, ~7 pcs/lb. 50092	\$13.25	
Size 10, ~5 pcs/lb. 50102	\$13.25	
Assortment, 1 lb. 50012	\$17.25	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$13.50

Non-fenestrated towels are sterile, latex-free and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515 \$25.50

Tape

Tape



Autoclave, steam


¾" x 60 Yd Roll. **44073ST048** \$8.50
 1" x 60 Yd Roll. **44073ST036** \$10.35↓

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll
 ½" Wide, each \$5.50
 Green. **134623005** Yellow. **134622005**
 ¾" Wide, Pk/4 \$23.75
 Green. **134820075** Yellow. **134810075**
 White. **134800075**
 3" Core, 40 yd roll
 ½" Wide, each \$10.05
 Green. **134633005** White. **134630005**
 ¾" Wide, each \$13.35
 Green. **134633075** Yellow. **134632075**
 ¾" Wide, White, Pk/4. **134830075** \$53.35
 ¾" Wide, Mixed, Pk/6. **134630600** \$81.20
 blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000 \$42.20



Fuses

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. 312-100	\$2.00
0.2A/250V. 312-200	\$2.00
1A/250V. 312-001	\$2.00
2A/250V. 312-002	\$2.00
3A/250V. 312-003	\$2.00
4A/250V. 312-004	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. 313-700	\$2.00
1A/250V. 313-001	\$2.00
2.5A/250V	\$2.00

5 x 20 mm

Fast-Acting, GMA


0.5A/250V. 235-500	\$2.00
1A/250V. 235-001	\$2.00
3A/250V. 235-003	\$2.00
5A/125V. 235-005	\$2.00

Time-Delay (Slow-Blow), Glass tube

0.5A/250V. 239-500	\$2.00
1A/250V. 239-001	\$2.00
2A/250V. 239-002	\$2.00
3A/250V. 239-003	\$2.00
5A/250V. 239-005	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. 215-001	\$2.00
-------------------------------	--------



Cone Point



Cup Point

Set Screws, steel

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. 77702	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. 77704	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. 77706	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. 77708	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. 77722	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. 77724	\$0.50

Tools




Cutter, diagonal, 5". PLR-766.00 \$10.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00 \$80.95

Gauges



Vernier Caliper measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

Vernier Caliper, digital, 3-button. 600VCD001 \$29.55



Compasses, replaceable tips, Germany

3" Long. DIV-320.00	\$11.95
6" Long. DIV-326.00	\$16.95
Tips (replacement), Pair. DIV-320.05	\$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A \$39.00



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

Tubing Sizer, metric diameter gauge, PP. 133900000 \$35.80

Continued