

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker. Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot®, permanent, extra fine point	ea	bx/12
Black. 44102	\$1.70	17.45
Blue. 44103	\$1.70	17.45
Green. 44105	\$1.70	17.45
Red. 44104	—	17.45



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Pens, Super ^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic		
Regular Tip (4 mm tapered nib). K-600	\$61.00	
(71310, 195505)		
Fine Tip (2.5 mm tapered nib). K-650	\$46.00	
(71312, 195506)		

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

China Marker (Wax Pencil), Bx/12. 02089	\$11.50
--	---------



Diamond PerfectMark™, permanent, black, Bx/10. 1000	\$23.20
--	---------

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

Fine Tip	ea	bx/12
Black. 36701	—	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	ea	bx/12
Black. 1735790	\$2.85	29.65
Blue. 1735792	—	29.65

Diamond-Tipped



Carbide-Tipped

Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$22.75
Diamond Tip (replacement). RDT-1	\$12.90




Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers

Small, 26 x 18 cm (10 x 7"). XH-2002	\$257.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$345.00

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

10-Slide, polystyrene, USA	ea	pk/3
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00



20-Slide, Simpport

Clear Lid (PETG), black base (ABS), each. M920-1	\$162.00
Black Lid (ABS), black base (ABS), each. M920-2	\$162.00

Continued