

## Stereo Microscopes



420T-430PLL-10

### National

**420-430PLL-10** Stereo, Binocular, 10-40X Zoom..... \$598.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 10X-40X total magnification; working distance of 79 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (LED with rheostat control and directional beam), bottom only (LED), or both top and bottom
- Post stand with locking support collar allows viewing specimens up to 93 mm tall
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Lifetime mechanical and optical/3-year electrical warranty

**420T-430PLL-10** Stereo, Trinocular, 10-40X Zoom..... \$692.00

- Trinocular version of 420-430PLL-10
- 2X SLR adapter & 0.4X C-mount adapter included

425-445PLL-10



### National

**425-445PLL-10** Stereo, Binocular, 7.5-50X Zoom..... \$959.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-50X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warranty


**425T-445PLL-10** Stereo, Trinocular, 7.5-50X Zoom..... \$1,023.00

- Trinocular version of 425-445PLL-10



Meiji ..... Access [www.meijitechno.com](http://www.meijitechno.com)  
and call us to place an order

## Stands



*Please call regarding availability of items not listed*

### Stereo Boom Stands, Diagnostic Instruments

*All boom stands require a focus mount*

- Standard, nylon bearings. **SMS16B** ..... \$599.00
- Ball Bearing
- Medium-Duty. **SMS6B** ..... \$759.00
- Heavy-Duty. **SMS20-\*** ..... *See note below*
- Articulating Arm. **SMS25** ..... \$1,409.00

*\*Access [www.spotimaging.com](http://www.spotimaging.com) and call us with microscope model and boom stand version to place an order*



Meiji PBH Stand  
with Adapter Sleeve

### Stereo Dual Illumination Stands, 6V/10W, Meiji

- Pole-Type Base, 84 mm focus mount. **PBH**..... \$563.00
- Arm-Type Base, 84 mm focus mount. **ABZH**..... \$612.00
- Some stereo brands require an adapter sleeve*
- Stereo Stand Sleeve for 76 mm dia. bodies ..... \$30.00
- Stereo Stage Plates, 94.5 mm
- Clear Glass. **MA569** ..... \$45.00
- Filter, Blue Frosted, 40 mm. **MA563/05**..... \$50.00
- Black & White Plastic. **MA568** ..... \$45.00