



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000**\$679.00

Stacking Platform. **BR2000-STACK**.....\$107.00

14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$38.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300**\$737.00

Stacking Platform. **BR1000-STACK**.....\$98.80

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$38.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**\$603.00

Double Platform, 3¼" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**\$647.00

Double Platform, 3¼" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100**\$318.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.



BT30

BT302

Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**\$716.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302**\$827.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3000

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3000**\$1,335.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**...\$229.00

250 mL Flask Magnetic Clamp. **H1000-MR-250**\$32.30

500 mL Flask Magnetic Clamp. **H1000-MR-500**\$32.30

1000 mL Flask Magnetic Clamp. **H1000-MR-1000**\$49.40



R4045

R4040



Tube Rotator, programmable, Benchmark (2-year warranty)

Roto-Bot Mini. **R4040**\$476.00

Roto-Bot. **R4045**.....\$728.00

- Speed adjusts from 1-80 rpm
- Tube capacity
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Time is programmable from 1 sec – 99 hrs, 59 min / Cont.