


HS120035/7



HS120044



Tube Storage Boxes, Polypropylene, Heathrow

- 81-Well (9 x 9), 13 x 13 x 4.7 cm
Holds 1.5-2.0 mL tubes and cryovials.
Features imprinted coordinates on a colored base and a translucent, friction-fit lid with frosted panels.
- Blue base, Pk/5. **HS120035** \$33.65
- Pink base, Pk/5. **HS120037** \$33.65
- Black base & lid, Pk/5. **HS120041** \$33.65
- 100-Well (10 x 10), mixed, Pk/5. **HS120044** \$35.10
Holds 1.5-2.0 mL tubes and cryovials.
Hinged lid has molded grid lines; base has imprinted and raised coordinates for each well. 14.1 x 15.1 x 5.7 cm.

Holds 1.5-2.0 mL tubes in a two-piece rack.




Tube Storage Boxes, Polystyrene Foam (styrofoam), USA

- 50-Well, Pk/10. **HS4283P** \$41.60
- 100-Well, Pk/10. **HS4282** \$52.05

HS2165C



HS2165A



HS2165B



HS2166




Tube Racks, Floating, HDPE Foam w/Handle, Heathrow

- 18-Hole, 0.2/0.5/1.5-2.0 mL tubes, blue circle, Pk/5. **HS2166** \$15.35
- 24-Hole, 1.5-2.0 mL tubes, yellow rectangle, Pk/5. **HS2165C** \$15.10
- 8-Hole, 15 mL tubes, blue parallelogram, Pk/5. **HS2165B** \$13.30
- 4-Hole, 50 mL tubes, green diamond, Pk/5. **HS2165A** \$14.85



Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

- 10-13 mm Tubes, 90-place, each. **187490000** \$26.25
- 15-16 mm Tubes, 60-place, each. **187490001** \$26.25
- 18-20 mm Tubes, 40-place, each. **187490002** \$26.25



Cube Rack




4-Way Rack




4-Way Rotatable Rack

Tube Racks, Interlocking, polypropylene, Heathrow

- Cube, mixed colors, Pk/5. **HS29050A** \$41.45
Holds (4) 50 mL tubes, (10) 15 mL tubes, (12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
- 4-Way, mixed colors, Pk/5. **HS29022G** \$38.20
Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes.
- 4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** \$40.30
Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.




HS24322A/E



HS120112

Tube Racks, Mega Rack, polypropylene, Heathrow

- 10-13 mm Tubes
216-Place, 33.7 x 26.7 x 9.1 cm(H)
Blue, each. **HS24322A** \$37.05
White, each. **HS24322E** \$37.05
- 432-Place, 67.4 x 26.7 x 9.1 cm(H)
Blue, each. **HS120113** \$68.00
White, each. **HS120112** \$68.00



80-Well



96-Well Reversible

Tube Racks, Microtube, polypropylene, Heathrow

- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
- Blue, Pk/5. **HS29025B** \$24.55
- Pink, Pk/5. **HS29025D** \$24.55
- Mixed, Pk/5. **HS29025G** \$24.55
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
- Blue, Pk/5. **HS2345C** \$47.85
- Green, Pk/5. **HS2345D** \$47.85
- Pink, Pk/5. **HS2345E** \$47.85

Continued