

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker. Pilot pens are ideal for cytologists and pathologists.



| Pens, Pilot®, permanent, extra fine point | ea | bx/12 |
|---|--------|-------|
| Black. 44102 | \$1.70 | 17.45 |
| Blue. 44103 | \$1.70 | 17.45 |
| Green. 44105 | \$1.70 | 17.45 |
| Red. 44104 | — | 17.45 |



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

| Pens, Super ^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic | ea | bx/12 |
|---|---------|-------|
| Regular Tip (4 mm tapered nib). K-600 | \$65.00 | |
| (71310, 195505) | | |
| Fine Tip (2.5 mm tapered nib). K-650 | \$50.00 | |
| (71312, 195506) | | |

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Wax Pencil

Markers, USA

| | |
|--|---------|
| Wax Pencil (China Marker), Bx/12. 02089 | \$11.50 |
|--|---------|



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

| | |
|--|---------|
| Diamond PerfectMark™, permanent, black, Bx/10. 1000 | \$24.20 |
|--|---------|

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

| Fine Tip | ea | bx/12 |
|-----------------------------|--------|-------|
| Black. 36701 | — | 29.65 |
| Blue. 36703 | \$2.85 | 29.65 |
| Red. 36702 | \$2.85 | 29.65 |
| Ultra-Fine Tip | ea | bx/12 |
| Black. 1735790 | \$2.85 | 29.65 |
| Blue. 1735792 | — | 29.65 |

Diamond-Tipped



Carbide-Tipped

Scribes, retractable

| | |
|--|---------|
| Carbide-Tipped (sharpenable). RSC | \$15.30 |
| Carbide Tip (replacement). RCT-1 | \$6.20 |
| Diamond-Tipped. RSD | \$22.75 |
| Diamond Tip (replacement). RDT-1 | \$12.90 |




XH-2002

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers

| | |
|---|----------|
| Small, 26 x 18 cm (10 x 7"). XH-2002 | \$257.00 |
| Large, 63 x 20 cm (25 x 8"). XH-2001 | \$345.00 |

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

| 10-Slide, polystyrene, USA | ea | pk/3 |
|----------------------------------|---------|--------|
| Clear. 240-9000-010 | \$54.95 | 154.00 |
| Black. 240-9020-Z10 | \$56.10 | 157.00 |



M920-1

20-Slide, Simport

| | |
|---|----------|
| Clear Lid (PETG), black base (ABS), each. M920-1 | \$166.00 |
| Black Lid (ABS), black base (ABS), each. M920-2 | \$166.00 |