


Staining Dishes

Glass (flint)

- 20-Slide, holds slides horizontally in metal rack
 Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)
 Set (dish, rack & lid). **SDS-20** \$26.60
 Dish & Lid. **SD-20** \$16.20
 Rack only, stainless steel. **SR-20** \$13.30
- 30-Slide, holds slides vertically in metal rack
 Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)
 Set (dish, rack & lid). **SDS-30** \$36.35
 Dish & Lid. **SD-30** \$18.05
 Rack only, stainless steel. **SR-30** \$22.25
- 60-Slide, holds slides vertically in metal rack
 Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)
 Set (dish, rack & lid). **SDS-60** \$50.85
 Dish & Lid. **SD-60** \$27.70
 Rack only, stainless steel. **SR-60** \$28.60



Plastic Set

Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

- Plastic (polyoxymethylene [POM]), holds 25 slides
 Set (dish, rack & lid). **0990-99** \$50.90



Staining Jars, hold slides (up to 76 x 26 mm) vertically

- Coplin, holds 10 slides back-to-back ea cs/12
- Glass w/Ground Glass Cover, USA. **673** \$11.45 128.00
- Glass Cover, Fits 673 Jar. **673-LID** \$4.30 —
- ea cs/12
- Glass w/Screw-On Phenolic Cap, **USA. 674** \$11.45 128.00
- Phenolic Cap, Fits 674 Jar, black. **674-LID** ...\$4.35 —
- ea cs/12
- PP w/Screw-On Cap. **600P90101** \$4.65 45.85
- Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is flat.
- ea cs/12
- Wide Top, Glass w/Square Lid, **USA. 675** \$15.35 171.00
- holds 16 slides back-to-back
- Glass Lid, Fits 675 Jar. **675-LID** \$4.65 —




Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Staining Jars (cont)

Wide Top, Plastic (POM)

- 12-Slide, EasyDip™, Simport ea cs/6
- Jar only, Blue. **M900-12B** \$10.80 59.30
- Jar only, Green. **M900-12G** \$10.80 59.30
- Jar only, Pink. **M900-12P** \$10.80 59.30
- Jar only, White. **M900-12W** \$10.80 59.30
- Jar only, Yellow. **M900-12Y** \$10.80 59.30
- Rack/Lid only, gray. **M905-12DGY** \$12.00 65.95



Anodized aluminum staining jar holder

(Jars sold separately)

- Staining Jar Holder, 6-Place. **M906** \$72.25



- 25-Slide, light green, microwaveable, acetal
- Jar & Lid. **4501030** \$11.55
- Rack only. **4501042** \$18.90
- Staining Jar Station, 12-Place. **4501011** \$221.00 ↓
- Includes metal frame, (12) 25-slide dishes w/ lids (1030) & (1) rack (1042).
- Dimensions: 29(L) x 5(W) x 5"(H)