



**Tissue (Embedding) & Biopsy Cassettes**



**Base Molds**

Plastic, disposable, Simport



- 15 x 15 x 5 mm, Cs/1000. **M475-2** ..... \$132.00
- 24 x 24 x 5 mm, Cs/1000. **M475-3** ..... \$132.00
- ★ 30 x 24 x 5 mm, Cs/1000. **M475-4** ..... \$132.00
- 37 x 24 x 5 mm, Cs/1000. **M475-5** ..... \$132.00
- 37 x 24 x 10 mm, Cs/500. **M475-10** ..... \$68.50

Designed for Simport Macrosette™.

Metal, reusable

- 15 x 15 x 5 mm, Pk/12. **HB-15** ..... \$23.00
- 24 x 24 x 5 mm, Pk/12. **HB-24** ..... \$23.00
- 30 x 24 x 5 mm, Pk/12. **HB-30** ..... \$23.00
- 37 x 24 x 5 mm, Pk/12. **HB-37** ..... \$23.00

Solvent resistant,  
0.2 x 0.2 mm mesh.




**Biopsy Bags**

- Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011** ..... \$235.00
- Paper (Tea Bags), 2¼ x 2¼", Pk/500. **56100** ..... \$40.55

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



**Biopsy Pads**


Rectangular

- Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014** ..... \$42.00
- Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1** ..... \$194.00

Circular, Simport

- Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2** ..... \$236.00


**Biopsy Cassettes**




Economy Biopsy Cassettes

**Biopsy Cassettes, Detached Lids (packed separately)**

Economy, White, small square pores, Pk/500. **BC-0109W** ..... \$42.90



Simport Histosette I



Simport Micromesh

**Biopsy Cassettes, Attached Lids, acetyl polymer**

Histosette™ I, 1 mm square pores, Pk/1500 ..... 1 pk pk/3  
\$129.00 363.00  
30° front labeling area.


- Blue. **M491-6**                      Pink. **M491-3**
- Green. **M491-4**                      White. **M491-2**

Unisette™, hinged lids, Pk/1500 ..... 1 pk pk/3  
\$138.00 392.00  
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

- Blue. **M506-6**                      Pink. **M506-3**
- Green. **M506-4**                      White. **M506-2**

Micromesh™, 0.38 mm pores, Pk/1000 ..... 1 pk pk/3  
\$216.00 602.00  
45° front labeling area.

- Blue. **M507-6**                      Pink. **M507-3**                      Yellow. **M507-5**
- Green. **M507-4**                      White. **M507-2**



Compatible with Primera Signature™ and Tissue-Tek SmartWrite™ cassette printers.

**Biopsy Cassettes, Taped**


Micromesh™, 0.38 mm pores, Pk/2000 ..... \$361.00  
45° front labeling area.

- Blue. **M407-6T**                      Pink. **M407-3T**                      Yellow. **M407-5T**
- Green. **M407-4T**                      White. **M407-2T**

Continued

## Tissue Cassettes

Economy  
Tissue  
Cassettes



### Tissue Cassettes, Detached Lids, slotted pores

Economy, Pk/500 ..... \$42.90  
 Blue. **EC-0109B**                      White. **EC-0109W**  
 Lilac. **EC-0109L**

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ II, Simport, Pk/1000 ..... \$88.10	254.00	
45° front labeling area.		
Blue. <b>M485-6</b>	Pink. <b>M485-3</b>	
Green. <b>M485-4</b>	White. <b>M485-2</b>	

Uni-Cassettes™, **USA**, Pk/500 ..... \$60.45

Aqua. <b>6059</b>	Lime. <b>6062</b>
Blue. <b>6054</b>	Peach. <b>6064</b>
Burnt Orange. <b>6056</b>	Pink. <b>6052</b>
Chartreuse. <b>6061</b>	Red. <b>6060</b>
Gold. <b>6058</b>	Rose. <b>6063</b>
Green. <b>6053</b>	Tan. <b>6055</b>
Lavender. <b>6065</b>	White. <b>6050</b>
Lilac. <b>6057</b>	Yellow. <b>6051</b>



Simport



Simport  
Macrosette

### Tissue Cassettes, Attached Lids, slotted pores

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ I, snap-apart lids, Pk/1500 ..... \$129.00	363.00	
30° front labeling area.		

Blue. <b>M490-6</b>	Pink. <b>M490-3</b>
Green. <b>M490-4</b>	White. <b>M490-2</b>

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ II, integral lids, Pk/1500 ..... \$172.00	494.00	
45° front labeling area.		

Blue. <b>M492-6</b>	Pink. <b>M492-3</b>
Green. <b>M492-4</b>	White. <b>M492-2</b>

	<u>1 pk</u>	<u>pk/3</u>
Uniset™, hinged lids, Pk/1500 ..... \$138.00	392.00	
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)		

Blue. <b>M505-6</b>	Pink. <b>M505-3</b>
Green. <b>M505-4</b>	White. <b>M505-2</b>

Macrosette™, 2X normal depth, white, Pk/750. **M512** ..... \$182.00  
 Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all  
available metal lids.

### Tissue Cassettes, without Lids, acetyl polymer, premium

Capsettes™, **USA**, Pk/500 ..... \$41.60

Aqua. <b>6039</b>	Lime. <b>6042</b>
Blue. <b>6034</b>	Peach. <b>6044</b>
Burnt Orange. <b>6036</b>	Pink. <b>6032</b>
Chartreuse. <b>6041</b>	Red. <b>6040</b>
Gold. <b>6038</b>	Rose. <b>6043</b>
Green. <b>6033</b>	Tan. <b>6035</b>
Lavender. <b>6045</b>	White. <b>6030</b>
Lilac. <b>6037</b>	Yellow. <b>6031</b>

Simport, 45° front labeling area

Green, Pk/1500. **M480-4** ..... \$96.55

## Cassette Storage

**Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342** ..... \$121.00  
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
in 4 drawers, each divided into 12 rows.  
Cabinets can be stacked up to 8 high.  
Dimensions: 49.5 x 38 x 21 cm.




**Tissue Cassette Cabinet, Steel, 4-drawer. 485400** ..... \$250.00



Simport

Store 1500 cassettes or 990 embedding  
rings in 6 stackable drawers.  
High-impact polystyrene resists  
temperatures from -80°C to +80°C.  
Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6** ..... \$157.00



LB-38



LB-30

### Brushes

Flat Bristles, paint style, camel hair  
 ½" Wide, 6¼" total length. **LB-36** ..... \$3.90  
 1" Wide, 7" total length. **LB-37** ..... \$6.15  
 2" Wide, 7¾" total length. **LB-38** ..... \$5.30  
 Round Bristles, camel hair, 8¾" long. **LB-30** ..... \$10.25