


Tissue Cassettes

Economy
Tissue
Cassettes



Tissue Cassettes, Detached Lids, slotted pores

Economy, Pk/500 \$42.90
 Blue. **EC-0109B** White. **EC-0109W**
 Lilac. **EC-0109L**

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ II, Simport, Pk/1000 \$88.10	254.00	
45° front labeling area.		
Blue. M485-6	Pink. M485-3	
Green. M485-4	White. M485-2	

Uni-Cassettes™, **USA**, Pk/500 \$60.45

Aqua. 6059	Lime. 6062
Blue. 6054	Peach. 6064
Burnt Orange. 6056	Pink. 6052
Chartreuse. 6061	Red. 6060
Gold. 6058	Rose. 6063
Green. 6053	Tan. 6055
Lavender. 6065	White. 6050
Lilac. 6057	Yellow. 6051



Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ I, snap-apart lids, Pk/1500 \$129.00	363.00	
30° front labeling area.		

Blue. M490-6	Pink. M490-3
Green. M490-4	White. M490-2

	<u>1 pk</u>	<u>pk/3</u>
Histosette™ II, integral lids, Pk/1500 \$172.00	494.00	
45° front labeling area.		

Blue. M492-6	Pink. M492-3
Green. M492-4	White. M492-2

	<u>1 pk</u>	<u>pk/3</u>
Uniset™, hinged lids, Pk/1500 \$138.00	392.00	
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)		

Blue. M505-6	Pink. M505-3
Green. M505-4	White. M505-2

Macrosette™, 2X normal depth, white, Pk/750. **M512** \$182.00
 Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all
available metal lids.

Tissue Cassettes, without Lids, acetyl polymer, premium

Capsettes™, **USA**, Pk/500 \$41.60

Aqua. 6039	Lime. 6042
Blue. 6034	Peach. 6044
Burnt Orange. 6036	Pink. 6032
Chartreuse. 6041	Red. 6040
Gold. 6038	Rose. 6043
Green. 6033	Tan. 6035
Lavender. 6045	White. 6030
Lilac. 6037	Yellow. 6031

Simport, 45° front labeling area

Green, Pk/1500. **M480-4** \$96.55

Cassette Storage

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342 \$121.00
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
in 4 drawers, each divided into 12 rows.
Cabinets can be stacked up to 8 high.
Dimensions: 49.5 x 38 x 21 cm.




Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$250.00



Simport

Store 1500 cassettes or 990 embedding
rings in 6 stackable drawers.
High-impact polystyrene resists
temperatures from -80°C to +80°C.
Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6 \$157.00



LB-38



LB-30

Brushes

Flat Bristles, paint style, camel hair

½" Wide, 6¼" total length. **LB-36** \$3.90

1" Wide, 7" total length. **LB-37** \$6.15

2" Wide, 7¾" total length. **LB-38** \$5.30

Round Bristles, camel hair, 8¾" long. **LB-30** \$10.25



Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. **246530002** \$69.25
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** \$96.15

Hinged
Polystyrene



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80



Brushes

Bench, horse hair, 10" long. **BRS-969.10** \$9.20
 Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".
 Brass Wire, 9 3/4" long. **LB-40** \$8.70
 Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".
 Flat Bristles, paint style
 Camel Hair
 1/2" Wide, 6 3/4" total length. **LB-36** \$3.90
 1" Wide, 6 1/2" total length. **LB-37** \$6.15
 2" Wide, 7 3/4" total length. **LB-38** \$5.30
 Horsehair, heavy-duty bristles, metal handle. **33761880** \$0.85
 Round Bristles, camel hair, 8 1/2" long. **LB-30** \$10.25
 Swivel-Head, black nylon bristles. **LB-14** \$4.20
 Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany
 Brass. **BRS-295.00** \$12.20
 Nylon. **BRS-293.00** \$9.70
 Steel. **BRS-292.00** \$12.80

Burners



Burner, Alcohol, 120 mL, glass body, metal cap. **97-5313** \$6.80
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** \$5.00



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



Burners, Bunsen, natural gas, Humboldt Mfg, **USA**

Standard, flame stabilizer. **H-5870** \$22.45
 High Temperature, Meeker-type. **H-5550** \$78.90

Continued