

Balances

Portable



Model	DWP-1001
Capacity (g)	1000
Readability (g)	0.1
Modes	g, oz, ct, gn
Pan	14.6 cm (5.75") stainless steel
Power	4 AA/AC adapter
Price	\$72.00

Warranty: 1 year

Toploading

Benchmark
W3200-1200



Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$1,081.00	\$874.00

Warranty: 2 years

Analytical

Benchmark



Model	★ W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,610.00	\$1,708.00

Warranty: 2 years

Baths



BSH1002

Aluminum
Blocks

Dry Baths, digital, Benchmark (2-year warranty)

1-Block Capacity. BSH1001	\$386.00
2-Block Capacity. BSH1002	\$467.00
4-Block Capacity. BSH1004	\$795.00

Blocks sold separately – please specify configuration

Blocks, aluminum

Holds (48) 0.2 mL Tubes. BSW02	\$98.00
Holds (24) 0.5 mL Tubes. BSW05	\$98.00
Holds (24) 1.5 mL Tubes. BSW1500	\$98.00
Holds (24) 1.5/2.0 mL Tubes. BSW1520	\$98.00
Holds (20) 12/13 mL Tubes. BSW13	\$98.00
Holds (12) 15 mL Tubes. BSW15	\$98.00
Holds (5) 50 mL Tubes. BSW50	\$98.00



40 x 1.5/2.0 mL



15 x 50 mL



41 x 15 mL

Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. B2000-4	\$799.00
8 L. B2000-8	\$869.00
12 L. B2000-12	\$1,036.00



Economy

Polypropylene
Beakers



Premium

Beakers

Glass (borosilicate), Griffin low form, single scale, economy

10 mL, no graduations, each. 6102022910	\$1.85
50 mL, each. 7-440050	\$3.10
100 mL, each. 7-440100	\$3.40
150 mL	
Pk/12. 61020229150	\$46.30
Cs/48. 61020229150	\$156.00
250 mL, each. 61020229250	\$3.50
400 mL, each. 61020229400	\$4.40
600 mL, each. 61020229600	\$5.15
800 mL	
Pk/6. 600BG800	\$62.80
Cs/24. 600BG800	\$231.00
1000 mL, each. 600BG1000	\$8.65
1000 mL, Heavy-Duty, each. 600BG1003	\$14.00
2000 mL, each. 600BG2000	\$15.85
4000 mL, each. 600BG4000	\$43.70

Continued

Beakers (cont)

Polypropylene

Economy, Triple-Spout

- 50 mL, 5 mL graduations, Cs/100. **3640** \$18.80
- 100 mL, 10 mL graduations, Cs/100. **80093** \$28.60
- 250 mL, 20 mL graduations, Cs/100. **3642** \$42.40
- 1000 mL, 50 mL graduations, Cs/100. **80098** \$96.25

Premium, Griffin low form, single scale

- 50 mL, 10 mL graduations, each. **262100000** \$2.95
- 100 mL, 25 mL graduations, each. **262110000** \$2.80
- 250 mL, 50 mL graduations, each. **262130000** \$4.95
- 600 mL, 100 mL graduations, each. **262150000** \$7.85
- 1000 mL, 200 mL graduations, each. **262160000** \$7.95
- 2000 mL, 500 mL graduations, each. **262170000** \$19.20
- 4000 mL, 1000 mL graduations, each. **262190000** \$39.70



Roll Bench Liner

Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liners, white

- Roll, 50' x 20", Labmat™. **246750000** \$66.65
- Sheets,
 - Bel-Art, 20 x 24", Pk/100. **246750020** \$232.00
 - Ahlstrom, 46 x 57 cm (18 x 22½"), Pk/50. **4002-4657** \$81.20

Bottles



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

- Bottle, Balsam, glass, 45 mL (empty), USA. 616** \$18.35



- Bottle Carrier, BOD, holds 8 bottles, polypropylene**
- Bx/2. **600BODCAR** \$36.85




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

- LDPE (low-density polyethylene), Yorker cap
 - 30 mL (1 oz), Bg/12. **116370001** \$18.20
 - 60 mL (2 oz), Bg/12. **116370002** \$21.15
 - 125 mL (4 oz), Bg/12. **116370004** \$35.25
 - 185 mL (6 oz), Bg/12. **116370006** \$35.25

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

- 20 mL, each. **18145** \$7.00

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

- 30 mL (1 oz), Bx/12. **3104053** \$14.30
- Bottle Dropper Assembly (replacement)
 - 20/400, 7 x 76 mm Pipet, Pk/12. **3104053D** \$8.00
- 60 mL (2 oz), Bx/12. **3104054** \$22.55
- Bottle Dropper Assembly (replacement)
 - 20/400, 7 x 89 mm Pipet, Pk/12. **12574054D** \$8.50

French Square



Bottles, Dropping, Clear Glass (flint)

- 15 mL (1/2 oz), Bx/12. **3104062** \$20.50 ↓
- 30 mL (1 oz), Bx/12. **3104063** \$25.65 ↓
- Bottle Dropper Assembly (replacement)
 - 24/400, 7 x 75 mm Pipet, Pk/12. **3104075** \$11.35

Continued