




**Boxes, Glass Disposal**, cardboard, polyethylene bag included  
 Bench Model, 8 x 8 x 10", Pk/6. **246530002** ..... \$69.25  
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** ..... \$96.15




**Boxes, Plastic**, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
<b>H-1</b>	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
<b>H-1/11B</b>	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
<b>H-54</b>	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
<b>H-54/2916G</b>	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
<b>H-28/334B2</b>	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
<b>H-10</b>	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
<b>H-21/512B</b>	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80



**Brushes**

Bench, horse hair, 10" long. **BRS-969.10** ..... \$9.20  
 Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".  
 Brass Wire, 9 3/4" long. **LB-40** ..... \$8.70  
 Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".  
 Flat Bristles, paint style  
 Camel Hair  
 1/2" Wide, 6 3/4" total length. **LB-36** ..... \$3.90  
 1" Wide, 6 1/2" total length. **LB-37** ..... \$6.15  
 2" Wide, 7 3/4" total length. **LB-38** ..... \$5.30  
 Horsehair, heavy-duty bristles, metal handle. **33761880** .... \$0.85  
 Round Bristles, camel hair, 8 1/2" long. **LB-30** ..... \$10.25  
 Swivel-Head, black nylon bristles. **LB-14** ..... \$4.20  
 Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".



Remove debris, rust, or corrosion in hard to reach places.

Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany  
 Brass. **BRS-295.00** ..... \$12.20  
 Nylon. **BRS-293.00** ..... \$9.70  
 Steel. **BRS-292.00** ..... \$12.80

**Burners**



**Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313** ..... \$6.80  
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** ..... \$5.00



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



**Burners, Bunsen**, natural gas, Humboldt Mfg, **USA**  
 Standard, flame stabilizer. **H-5870** ..... \$22.45  
 High Temperature, Meeker-type. **H-5550** ..... \$78.90

*Continued*