



Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. **246530002** \$69.25
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** \$96.15

Hinged
Polystyrene



Boxes, Plastic, (see following chart)

| Plastic Boxes | | | | |
|-------------------|------------------------|-------------|---------------------------|--------|
| Catalog # | Inside Dimensions | Foam Insert | Suggested Application | Price |
| H-1 | 1 x 1 x 1/4" | - | Reticles, Pointers | \$0.40 |
| H-1/11B | 1 x 1 x 1/4" | 1 | Reticles, Pointers | \$0.50 |
| H-54 | 2 9/16 x 2 9/16 x 1/2" | - | Filters | \$1.20 |
| H-54/2916G | 2 9/16 x 2 9/16 x 1/2" | 1 | Filters | \$1.35 |
| H-28/334B2 | 3 3/4 x 1 3/16 x 5/8" | 2 | Slides, Stage Micrometers | \$1.40 |
| H-10 | 2 7/8 x 1 3/16 x 3/8" | - | Small Parts | \$0.60 |
| H-21/512B | 5 1/2 x 1 3/4 x 1/2" | 1 | Forceps & Scissors | \$1.80 |



Brushes

Bench, horse hair, 10" long. **BRS-969.10** \$9.20
 Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".
 Brass Wire, 9 3/4" long. **LB-40** \$8.70
 Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".
 Flat Bristles, paint style
 Camel Hair
 1/2" Wide, 6 3/4" total length. **LB-36** \$3.90
 1" Wide, 6 1/2" total length. **LB-37** \$6.15
 2" Wide, 7 3/4" total length. **LB-38** \$5.30
 Horsehair, heavy-duty bristles, metal handle. **33761880** \$0.85
 Round Bristles, camel hair, 8 1/2" long. **LB-30** \$10.25
 Swivel-Head, black nylon bristles. **LB-14** \$4.20
 Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany
 Brass. **BRS-295.00** \$12.20
 Nylon. **BRS-293.00** \$9.70
 Steel. **BRS-292.00** \$12.80

Burners



Burner, Alcohol, 120 mL, glass body, metal cap. **97-5313** \$6.80
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** \$5.00



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



Burners, Bunsen, natural gas, Humboldt Mfg, **USA**

Standard, flame stabilizer. **H-5870** \$22.45
 High Temperature, Meeker-type. **H-5550** \$78.90

Continued



Burner Lighter, replaceable single flint, round file. **SOL-290.00** ... \$3.40
 Flint (replacement), Pk/5. **SOL-290.01** \$1.85

Carboys & Jugs



Carboys with Spigot, polypropylene, round
 10 L (2.5 gal). **118460025** \$70.20
 20 L (5 gal). **60034052** \$97.95



Jug, dispensing, HDPE jug/LDPE spigot
 20 L (5 gal). **118500000** \$46.75



Carts, utility, stainless steel, Gusto
 300 lb Capacity. **CART300S** \$273.00
 26.75"(L) x 15.5"(W) x 36"(H)
 500 lb Capacity. **CART500S** \$409.00
 29.75"(L) x 18"(W) x 36"(H)



Centrifuges, Mini
 5,500 RPM (2000 x g), 12-place, MyFuge™ 12. **C1012** \$334.00
 (1-year warranty)
 6,000 RPM (2000 x g), 8-place, MyFuge™ Mini. **C1008-*** \$246.00
 *Specify clear (-C), green (-G), or purple (-P) for lid color.
 (1-year warranty)
 10,000 RPM (5000 x g), 6-place, blue body. **MLX-110** \$193.00
 (2-year warranty)

Clamps & Supports

Rectangular
Steel Base




Triangular
Cast-Iron Base

Support (Ring) Stands, finished metal base/plated support rod

Rectangular Steel Base
 6 x 4" base, 5/16 x 18" rod. **3104SSS66** \$13.15
 8 x 5" base, 3/8 x 20" rod. **3104SSS67** \$16.10
 9 x 6" base, 1/2 x 24" rod. **3104SSS68** \$18.90

Rectangular Cast-Iron Base
 Provides greater stability than steel base.
 6 x 4" base, 5/16 x 18" rod. **3104SSC14** \$14.60
 8 x 5" base, 3/8 x 20" rod. **3104SSC15** \$17.20
 9 x 6" base, 1/2 x 24" rod. **3104SSC16** \$23.40

Triangular Cast-Iron Base
 5" Legs, 1/2 x 24" rod. **1257G2324** \$25.00
 6" Legs, 1/2 x 36" rod. **1257G23** \$30.65



Support Clamps
 Burette, plastic coated round jaws, 360° rotatable. **3104BC37** \$6.90
 Burette, Double, deluxe. **600CLBURDD** \$18.80
 Plastic covered spring jaws on a metal frame securely grip
 burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.
 Hook Collar, nickel plated. **1257G49** \$5.60
 Fits rods up to 1/2" (13 mm) in diameter.
 Right-Angle, aluminum. **97-4023** \$6.10
 Fits rods up to 5/8" in diameter.
 Thermometer, aluminum, rotates 360°. **600THCL01** \$7.65
 Fits rods up to 1/2" (13 mm) in diameter.



Support Rings, plated cast-iron
 3" OD (2 1/2" ID). **600SRC103** \$5.90
 4" OD (3 1/2" ID). **600SRC104** \$7.70
 5" OD (4 1/2" ID). **600SRC105** \$8.15

Continued